
EasyDrill v2[®]

Cranial Perforator

Broca para Perfuração Cranial EasyDrill V2
EasyDrill V2 Cranial Perforator
EasyDrill V2 Broca para Perforación Craneal

M I C R O M A R

www.micromar.com

CONTENTS

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INTENDED USE: The EasyDrill V2 Cranial Perforator is a sterile, single-use cutting device intended for performing cranial bone trephination.

CONTRAINDICATIONS: The perforator must not be used on bone thicknesses smaller than those described in the table located on the last page of this manual. Failure to comply could result in serious harm to the patient.

DEVICE DESCRIPTION: The EasyDrill V2 Cranial Perforator is a bone-drilling and -cutting device used in conjunction with a surgical motor and with a Hudson-type coupling attached to a speed-reduction device, in order to drill holes for access to the skull of patients of all ages whose treatment and/or diagnosis requires access to a certain zone of the brain. When properly used, the EasyDrill V2 Cranial Perforator opens a hole to allow a number of neurosurgical procedures through a decoupling mechanism, automatically stopping after the drilling is done and as soon as the drill no longer feels the resistance caused by the bone. Drilling the hole results in minimal exposure, which reduces trauma not only for the tissues around the brain, but also – and most importantly – for the brain itself.

WARNINGS

- Failure to follow this Instructions for Use manual and of the associated use equipment can lead to incorrect use of the device and therefore misuse, which increases the risk of death to the patient.
- To avoid involuntary movements during the surgical procedure, it is recommended that the patient's head be properly secured to a skull fixation device, which in turn should be properly fixed to the operating table.
- Do not use the drill in a hole or an area that has already been drilled, or the surrounding areas. The composition of reconstructed bone differs from that of healthy bone, increasing the risk of unwanted breakage. Leave a minimum spacing of 15 mm from one hole to another.

Note: The device may stop and/or shut down in the middle of trephination, there is the likelihood that the device may not reengage.

- It is essential to keep the drill perpendicular (90°) (according to Figure 3) at the predetermined point of the skull to be drilled, as an excessive deviation from perpendicularity may cause the product to fail and lead to serious patient injury.
- Select a drill bit suitable for the thickness according to the table on the last page. This prevents the drill from tearing the bone or brain tissue (similar effect when the bit is not positioned at 90°).
- If the components of the EasyDrill V2 Cranial Perforator become loose during trephination, discontinue use and contact Micromar or its representative.
- Take extra care with the pressure applied when drilling the skull of infants, children, and the elderly or pathological bone.
- The drill can cut or tear the dura mater when it unlocks. Check some conditions before performing the procedure, such as the existence of adhering dura mater or other anomalies adjacent to the drilling site.

PRECAUTIONS

- Do not use the EasyDrill V2 Cranial Perforator if the sterile packaging is opened or damaged.

- Do not use the EasyDrill V2 Cranial Perforator if the label is damaged or missing from the outer plastic drive body.
- The EasyDrill V2 Cranial Perforator is a single-use cutting device, which should never be used on more than one patient. Biologically hazardous drilling and cutting instruments must be safely discarded according to the applicable local regulations.
- This product must be handled by a surgeon specialized in neurosurgery.
- The surgeon must assess the patient's condition, including his/her age, weight, and health status before planning the surgery.
- Do not re-use, reprocess, or re-sterilize this product. Using an EasyDrill V2 Cranial Perforator device that has already been used may cause injury to the patient, operator and/or operating room staff; Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

PRIOR TO USE: SYSTEM SET UP AND FUNCTIONAL TEST

System Setup

- The EasyDrill V2 Cranial Perforator should only be used with a motorized driver with a Hudson-type chuck.
- The EasyDrill V2 Cranial Perforator must be used with the drilling equipment's motor rotation in clockwise direction.

Note: Verify your drilling equipment's type of connection and request an adapter if it is not of the Hudson-type.

EasyDrill V2 Cranial Perforator Functional test

Perform the following function test before each trephination.

Use sterile disposable gloves to perform the test.

CAUTION: If the sterile protective cover is not available, use sterile gauze (as a protective barrier) when handling the cutting end of the EasyDrill V2 Cranial Perforator.

NOTE: Do not use the product and contact customer service if either of the two function tests below fail to present the indicated results:

1. Firmly hold the EasyDrill V2 Cranial Perforator's plastic drive shank (Hudson shank) and the protective cover, turning the cutting end of the EasyDrill V2 Cranial Perforator device in a clockwise direction. The cutting end of the device should turn freely; see Figure 1.
2. Next, grasp the EasyDrill V2 Cranial Perforator by the plastic drive shank (Hudson shank) and press the protected cutting end of the device into the opposite hand, turn the drill clockwise until it locks; see Figure 2; to unlock, cease to apply pressure.

CRANIAL BONE DRILLING PROCEDURE:

Note: Recommended speed range for drilling using the EasyDrill V2 Cranial Perforator is 800 to 1000 RPM.

- After the functional tests, attach the EasyDrill V2 Cranial Perforator to the Hudson-type chuck of the motorized driver. Verify that the EasyDrill V2 Cranial Perforator is firmly fixed to the Hudson-type chuck before proceeding.
- Remove the protective cover from the cutting end of the EasyDrill V2 Cranial

Perforator device before drilling.

- Position and press the cutting end of the EasyDrill V2 Cranial Perforator gently against the skull to be perforated.
- Manually turn the drill until you feel the engagement movement.
- Position the EasyDrill V2 Cranial Perforator perpendicular to the skull (90°) on the point of perforation (see Figure 3). Supply power to the drilling system and begin trepanation.
- Maintain a gentle and continuous pressure while drilling.
- Irrigate the point of perforation with physiological serum while drilling to prevent bone necrosis.
- At the drilling depth where the EasyDrill V2 Cranial Perforator device ceases to find resistance to bone, it automatically disengages the cutting end of the device which stops rotating.

Note: At this stage, the plastic drive shank of the device will appear to continue to rotate.

Note: At this point, discontinue powering the system.

- Carefully remove the EasyDrill V2 Cranial Perforator from the skull.
- There may still be a thin bone layer in the hole that must be carefully removed with appropriate surgical instrument. If this layer is adhered to the brain material, take appropriate care not to injure the dura.
- Before performing another perforation, wipe the cutting end of the EasyDrill V2 Cranial Perforator removing the residues from the previous drilling. Perform testing procedures again.

LIMITED WARRANTY:

Micromar Ind. Com. Ltd. warrants this device to be free from design, and manufacturing defects. The device has a guarantee, provided by Micromar Ind. Com. Ltd. provided the following conditions are met:

- 1.1. The device is used in accordance with the Instructions for Use;
 - 1.2. The device is single used, in one patient, in the same surgical procedure.
- The limited warranty will be excluded for any of the following causes:
- 1.3. The device has been used inappropriately, or was subjected to misuse, or the Instructions for Use were not adequately followed;
 - 1.4. The device has been reused or reesterilized;
 - 1.5. The device has been used after the expiration date.

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