

# Generation of anti- pHLA antibodies

Subtitle The Key to High

*First name Last name, Senior Research Scientist*

Version no., Date etc



TCR-MIMETIC DISCOVERY · PHLA TARGETING

## Discriminating mutant from wild-type pHLA at the antibody discovery stage

**4**

 NAÏVE HUMAN LIBRARIES  
SCREENED IN PARALLEL

**53**

 MUTANT-PHLA-SPECIFIC  
UNIQUE FVS RECOVERED

**82%**

 MUTANT-SPECIFIC  
CLONES  
(BEST LIBRARY)

**2.97 nM**

 LEAD TCE AFFINITY  
KD ON MUTANT PHLA

### BACKGROUND


#### Targeting intracellular oncoproteins via the HLA class I axis

The therapeutically relevant proteome is dominated by intracellular targets that lie beyond the reach of conventional antibodies. HLA class I mediated antigen processing, however, surfaces fragments of these proteins as peptide–HLA (pHLA) complexes, generating substrates that are accessible to both T-cell receptors (TCRs) and TCR-mimetic (TCRm) antibodies. Mutation-driven oncoproteins are particularly attractive targets within this paradigm: their wild-type (WT) counterparts are presented across healthy tissues, often differing from the mutant peptide by a single residue. A TCRm therapeutic must therefore discriminate between sequences whose only difference may be one amino acid – failure to do so converts on-target activity into systemic on-target/off-tumour toxicity. The therapeutic index for these modalities is set at the discovery stage, not later in development.

### STRATEGY

#### Stringent counter-selection against WT and irrelevant pHLA

Phage display selections were performed against the mutant pHLA target across four FairJourney Bio fully human naïve libraries (each  $>1 \times 10^{10}$  functional diversity; scFv and Fab formats; separated kappa and lambda repertoires). Counter-selection pressure was applied during each round at increasing concentrations of irrelevant and WT pHLA, depleting cross-reactive phage before positive selection on the mutant target (Figure 1). Stringency was titrated empirically: as counter-selection concentration rose, fold-enrichment dropped sharply for libraries skewed toward non-specific binders, while libraries enriched in genuinely mutant-discriminating clones retained signal – providing an early read on which library outputs were worth carrying into screening.



**Figure 1.** End-to-end workflow. Four naïve human libraries (Fab and scFv) feed phage display selections under graded counter-selection pressure with irrelevant and WT pHLA. Output is triaged by ELISA-based specificity screening, HCDR3 sequence analysis, and FACS confirmation on T2 cells. The lead Fv is reformatted as a T-cell engager (TCE) and fully characterised in-house – affinity by SPR, function in cell-based assays.


### RESULTS

#### Discovery Output – 53 Mutant-Phla-Specific Fvs Across Four Libraries

All four libraries returned mutant-pHLA-specific clones, but with markedly different specificity yields and uniqueness profiles (Table 1). Library 3 emerged as the cleanest output: 82% of screened clones bound the mutant pHLA in ELISA with negligible irrelevant-pHLA cross-reactivity, resolving to 17 unique Fvs and 16 unique HCDR3 sequences. Library 4, run at lower counter-selection stringency, returned 71% mutant-specific but with 26% irrelevant-pHLA binding still present – illustrating directly the cost of insufficient counter-selection pressure. FACS on T2 cells pulsed with mutant peptide, irrelevant peptide, or unpulsed confirmed cellular-level specificity for the lead clones (Figure 2). 53 unique mutant-pHLA-specific Fvs were taken forward.

Library	Clones screened	Mutant pHLA	Irrelevant pHLA	Specific clones	Unique Fv	Unique HCDR3
Library 1	108	65%	6%	70	11	10
Library 2	108	70%	0%	76	8	3
Library 3	122	82%	7%	100	17	16
Library 4	122	71%	26%	87	17	15
<b>Total</b>	<b>460</b>	<b>—</b>	<b>—</b>	<b>333</b>	<b>53</b>	<b>44</b>


**Table 1.** Per-library specificity yield from ELISA screening on mutant vs. irrelevant pHLA, with downstream sequence diversity. Library 3 (highlighted) returned the highest specificity yield and the largest pool of unique HCDR3 sequences.



**Table 2.** Hit triage. **Left:** per-library hit-rates by ELISA on mutant vs. irrelevant pHLA. **Centre:** HCDR3 diversity and ELISA OD450 nm distribution per unique sequence. **Right:** FACS on T2 cells pulsed with mutant peptide (PC), irrelevant peptide (NC), or unpulsed – confirming cellular-level specificity.

### Lead TCE – Nanomolar Affinity, Mutant-Restricted T-Cell Killing

The lead Fv was reformatted as a T-cell engager (anti-pHLA × anti-CD3 architecture) and characterised end-to-end in-house. Affinity to the mutant pHLA was measured by SPR single-cycle kinetics across a 6–200 nM range, returning  $K_D = 2.97$  nM under a 1:1 Langmuir fit. The TCE drove dose-dependent IFN- $\gamma$  release and target-cell killing against HLA-A2+ cells presenting the mutant target at a 2:1 effector-to-target ratio. Specificity was tested across a six-cell-line panel – three lines presenting the WT target, three presenting the mutant target. Cytotoxicity and IFN- $\gamma$  release were strictly confined to the mutant-bearing lines: WT-presenting lines remained at baseline IFN- $\gamma$  across the full TCE dose range tested (1–10 nM), confirming that the mutant-versus-WT discrimination established at the Fv selection stage was preserved in the final TCE format (Figure 3).



**Figure 3.** TCE characterisation. Left: Single-cycle SPR sensorgram of TCE binding immobilized mutant pHLA, with the 1:1 Langmuir fit overlaid. Fivestep injection ladder spanning the 6.25–200 nM range;  $K_D = 2.97$  nM. Centre: dose-dependent IFN- $\gamma$  release (A) and cytotoxicity (B) against HLA2+ mutant-target-bearing cells (cell line 6) at 2:1 effector-to-target ratio. Right: TCE activity confined to mutant-target-bearing lines (cell lines 4–6); no specific activation on WT-target-bearing lines (cell lines 1–3) at TCE doses up to 10 nM. Significance: \*\*\*\*  $p < 0.0001$ .

## DISCUSSION

---

### Implications for safe TCRm therapeutic development

For mutant-pHLA targets, the therapeutic window is set at the moment of selection: any cross-reactive clone admitted into the discovery output will, in a TCE format, drive cytotoxicity in healthy tissue. Counter-selection at the phage-library stage is the cleanest place to impose that filter, because off-target binders are removed before commitment to a candidate is made. Two further considerations follow. First, peptide-level cross-reactivity to sequence-related but non-target pHLAs cannot be excluded by WT counter-selection alone – the short length of MHC-I peptides means that residues outside the mutation hotspot may permit cross-recognition of unrelated self-peptides. Downstream off-target panels (broader healthy-issue pHLA repertoires; in-silico cross-reactivity prediction) remain essential for clinical stage candidates. Second, format flexibility matters: the same Fv must retain its specificity profile across IgG, bispecific, and TCE architectures. The end-to-end workflow demonstrated here – naïve-library discovery, counter-selection-driven specificity, and in-house TCE reformatting and characterisation – closes that gap by maintaining a single experimental thread from library to functional candidate.

#### BOTTOM LINE

Discovery-stage counter-selection against WT and irrelevant pHLA, executed across four fully human naïve libraries and confirmed by orthogonal screening, delivers TCRm Fvs that retain mutant-versus-WT discrimination through TCE eformatting and into functional cell-killing assays – the prerequisite specificity profile for translating a TCRm modality into the clinic.

---

#### Discuss your TCR-mimetic discovery programme with our team:

[info@fjbio.com](mailto:info@fjbio.com) | [fjbio.com](http://fjbio.com)

Porto · Cambridge · San Francisco · San Diego



[www.fjbio.com](http://www.fjbio.com)  
[info@fjbio.com](mailto:info@fjbio.com)