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
ESTEREOTÁXICO AIMSISTEM TM03B

- ACESSÓRIOS

M I C R O M A R

MANUAL DO USUÁRIO	3
USER MANUAL	39
MANUAL DEL USUARIO	75

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ESTEREOTÁXICO AIMS SYSTEM TM03B

ACESSÓRIOS

ÍNDICE

1. INTRODUÇÃO 6

2. ILUSTRAÇÕES 6

3. COMPATIBILIDADE 6

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5. SEGMENTO "A" - TM0021FAA 7

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4. USE OF ACCESSORIES.....42

5. "A" SEGMENT - TM0021FAA.....43

6. LOCALIZERS FOR NUCLEAR MAGNETIC RESONANCE RM1000FAA / RM0002FAA.....44

7. EAR ALIGNERS - TM0004FAA.....45

8. RETRACTABLE ARMS SUPPORT - OP1010FAA.....51

9. SUPPORTS FOR COMPUTERIZED TOMOGRAPHY.....53

10. MANUAL CALCULATION FOR STEREOTACTIC COORDINATES.....58

11. TRAYS FOR TRANSPORT / STORAGE / STERILIZATION.....66

12. CLEANING AND STERILIZATION.....66

13. LUBRICATION.....68

14. PREVENTIVE AND CORRECTIVE MAINTENANCE.....69

15. SYMBOLS.....70

16. CARE AND PRECAUTIONS.....71

PRODUCT DESCRIPTION:

STEREOTACTIC AIMS SYSTEM TM03B

MODEL: TM03B - ACCESSORIES

ONLY USE THE EQUIPMENT AFTER CAREFULLY READING THE USER'S MANUAL AND RESPECTIVE WARNINGS.
CHECK IF THIS VERSION INSTRUCTION FOR USE IS ACCORDING TO THE PRODUCT PURCHASED.

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INFO@MICROMAR.COM

WWW.MICROMAR.COM

MICROMAR IND. COM. THE RIGHT TO MAKE CHANGES TO THIS MANUAL WITHOUT PRIOR NOTICE IS HEREBY RESERVED. CONTACT OUR SALES DEPARTMENT TO RECEIVE THE LATEST VERSION OF THE USER MANUAL.

THIS MANUAL CONTAINS INFORMATION PROTECTED BY COPYRIGHT

THIS MANUAL MAY NOT BE REPRODUCED, TRANSLATED OR TRANSMITTED BY ANY ELECTRONIC OR OTHER MEANS, IN WHOLE OR IN PART, WITHOUT EXPRESS AUTHORIZATION FROM MICROMAR

1. INTRODUCTION

This manual aims to provide assist the user on the safe and effective use the Stereotactic AimSystem TM03B.

Before using the equipment (TM03B and its accessories), the user should carefully read the manual and gain a total understanding of its contents, in addition to undergoing the training provided by Micromar. This manual should be kept with the equipment for reference in the event that questions arise.

2. ILLUSTRATIONS

The images in this manual show an illustration of the equipment and accessories. There may be some differences compared to prior equipment models and versions. These differences in illustrations do not interfere in its use.

3. COMPATIBILITY

The Accessories described in this manual may only be used in combination with other equipment or components which are expressly recognized by Micromar as compatible. Contact Micromar for more information..

4. USE OF ACCESSORIES

The Stereotactic AimSystem TM03B and its Accessories are made to be assembled and used according solely to the procedures indicated in this manual and for the purpose for which it was developed.

5. "A" SEGMENT - TM0021FAA

- **Recommended use**

The "A" Segment is used in conjunction with the Posterior Extension Posts and is recommended when indirect fixation needs to be done using three screws (two posterior and one anterior), allowing for approaches on the frontoparietal region.


Figure 1

- **"A" Segment Assembly**

1. To assemble the "A" Segment on the Frame, the Straight or Curved Segment must first be removed;
2. After removing the Straight Segment or Curved Segment, insert the "A" Segment so that it easily fits;
3. After the "A" Segment has been placed, tighten the screws using the 4mm Allen Key, turning them clockwise. The screws are placed on the bottom part of the Frame;
4. To remove the "A" Segment, loosen the screws securing it by using the 4mm Allen Key, turning it counterclockwise


Figure 2

- ⚠ **WARNING:** Fixation of "A" Segment screws should be done on the bottom part of the Frame;
- ⚠ **WARNING:** Never change the "A" Segment when the Frame is attached to the patient's head, since this will cause irreversible damage the Frame;
- ⚠ **WARNING:** When the Stainless Steel Screws are used, the Thread Adapter Bushing Guides need to be used.

6. LOCALIZERS FOR NUCLEAR MAGNETIC RESONANCE RM1000FAA / RM0002FAA

• Recommended use

The Localizers for Magnetic Resonance help the user to obtain images using the Nuclear Magnetic Resonance exam.

⚠ WARNING: Localizers for Magnetic Resonance can only be used with the Curved or Straight Segment;;

⚠ WARNING: When the Curved Segment is used, the Frame should be attached backwards, so that the Curved Segment is downwards;;

⚠ WARNING: Localizers for Magnetic Resonance can only be used in conjunction with the Extension Posts;

⚠ WARNING: When the Extension Posts are used, the localizers should always be in the "UPWARD" position;

⚠ WARNING: Incorrect placement of localizers will result in calculation errors in the neurosurgical procedure. Closely follow the instructions for use of Localizers for Magnetic Resonance;

• Assembly of Localizers for Magnetic Resonance

1. Open the velcro that attaches the NMR Localizers on the package;
2. Remove the localizers and Back Support from the package;
3. Fill the holes in the localizers with contrast liquids such as Gadolinium or Boric Water, ensuring that there are no bubbles in the holes' openings;
4. Assemble the localizers in their respective positions and the Back Support onto the Frame, locking them into place by using the localizer set screws;
5. When the Extension Posts are used with the Curved Segment, ensure that the Frame is assembled with the Curved Segment facing down;
6. Perform the imaging exam.


Figure 3

⚠ WARNING: Calculations for the resonance exam use the same method as computerized tomography.

7. EAR ALIGNERS - TM0004FAA

- **Recommended use**

Ear aligners are used to help in symmetrical fixation of the Frame to the patient's head.

- **Assembly of ear aligners using the straight segment**

1. Insert the body of the housing (right and left) into the Rulers;


Figure 4

2. Position the housing in the desired position;


Figure 5

3. Turn the fixation screws on the body of the housing clockwise to prevent movement;


Figure 6

4. Insert the Frame set;


Figure 7

5. Use the fixation screws found in the CT support to lock the aligner set into place, turning them clockwise;


Figure 8

6. Loosen the screw on the side of the housing and insert the ear aligner into the housing, adapting it according to the patient's ear anatomy;


Figure 9

7. Insert the Frame over the patient's head, positioning it in the region of interest;


Figure 10

8. If necessary, change the coordinates on the Ruler (on both sides) to make the Frame symmetrical and lock it in to place by turning the fixation screw clockwise;


Figure 11

9. Insert the Frame fixation screws;
10. Tighten the Frame fixation screws;
11. Remove the ear aligners.

- **Assembly of the Ear Aligners by using the Extension Posts and the Curved Segment or “A” Segment”**

1. Insert the body of the housing (right and left) into the Millimetric Rulers;


Figure 12

2. Position the housing in the desired position;


Figure 13

3. Turn the fixation screws on the body of the housing clockwise to prevent movement;


Figure 14

4. Insert the aligner set into the Frame;


Figure 15

5. Use the Fixation Screws found in the CT Support to lock the aligner set into place, turning them clockwise;


Figure 16

6. Loosen the screw on the side of the housing and insert the Ear Aligner into the housing, adapting it according to the patient's ear anatomy;


Figure 17

7. Insert the Frame over the patient's head, positioning it in the region of interest;


Figure 18

8. Lock the ear aligners in place by using the fixation screw on the side of the housing;


Figure 19

9. Insert the Frame fixation screws;
10. Tighten the Frame fixation screws;
11. Remove the Ear Aligners.

8. RETRACTABLE ARMS SUPPORT - OP1010FAA

- **Recommended use**

The Retractable Arms Support is a mechanical device that helps to adapt the Micromar Stereotactic system's retractable arms.

- **Assembly of Retractable Arms Support**

1. Assemble the Micromar Stereotactic System using the desired coordinates;
2. Check that the fixation screws are loose;


Figure 20

3. Insert the Support into the Frame


Figure 21

4. Use the 6mm Allen Key to lock the Support into place on the Frame, turning it clockwise;


Figure 22

5. Release the Handle that locks the retractor arm, turning it counterclockwise;


Figure 23

6. Insert the arm into the connection profile;


Figure 24

7. Turn the Handle clockwise to lock the arm in place.


Figure 25

NOTE: This support is only recommended for Leyla-type Retractable Arms.

9. SUPPORTS FOR COMPUTERIZED TOMOGRAPHY


Figure 26

- Recommended use**

The Supports for Computerized Tomography help the user to obtain images using the tomography exam.

Item	Model
A	UNIVERSAL (GE8800, SHIMADZU, SIEMENS, ELSINT, TOSHIBA, PICKER, PACE)
B	ELSCINT EXCEL 2400, 820
C	GE 8800, SRI, PROSPEED 9800, CTMAX 640, CTMAX OLD
D	GE, PACE, SYTEC 1800, 3000, 400, SYNERGY PLUS
E	GE 9800
F	GE HISPEED LXI, SX POWER, LIGHT SPEED ULTRA, LIGHT SPEED 8, 16, 32 SLICES
G	HITACHI W450, W4-40, W500, W2000, XR6000, PHILLIPS CX-S, CX-Q, SR4000 E AV
H	SHIMADZU (ALL)
I	SIEMENS SOMATOM DR, PHILLIPS MX800, BRILLIANCE CT 6, 10, 16 (UNTIL 10/2004)
J	PHILLIPS ALL MODELS UNSPECIFIED TOMOSCAN
K	PHILLIPS CT16,40 BRILLIANCE (UNTIL 11/2004)
L	PHILLIPS AURA
M	TOSHIBA TCT300S, TCI300S, TCI500S, 600HQ, 600XT
N	TOSHIBA AQUILION/ASTEION
O	PICKER PQ 2000, PQ 5000

⚠ WARNING: *There is a support model for each type of tomography.*

- Assembly of support for CT**

1. Remove the Handle that attaches the CT localizers, turning it counterclockwise;

Figure 27

2. Remove the CT Localizers;


Figure 28

3. Insert the support for CT on the tomography table;
4. Attach the support for CT (according to the support model or table configuration);
5. Release the blocking lever;


Figure 29

6. Release the Handles that attach the Frame to the support, turning them counterclockwise;


Figure 30

7. Insert the Frame into the support for CT;
8. Place the Handles in the Downward-Upward or Upward-Downward position;


Figure 31

NOTE: If the target is located above the Frame (Upward), the Handles should be inserted in the Upward position, and if the target is located below the Frame (Downward), the Handles should be inserted in the Downward position

9. Turn the Handles clockwise to attach the Frame;
10. Adjust the angular position of the support;


Figure 32

11. Lock the blocking lever.

- **Assembly of localizers for computerized tomography**

12. Place the Localizers in the direction of the target (Upward or Downward);
13. Use the fixation screws to lock the Frame Localizers into place;


Figure 33

14. In cases where the standard localizer cannot be used (obese patients or children), use the Short Localizer, ensuring that the Localizer's length will suit the intended field of vision;


Figure 34

15. The Radiological Bulkheads can also be used with the Extension Posts and "A" and Curved Segments.

WARNING: To correctly use the lateral (right and left) localizers, they must always be assembled in such a way

as to form a diagonal "N_i;"


Figura 35

- ⚠ WARNING:** When using the Extension Posts, the radiological bulkheads will always be in the "Upward" position;
- ⚠ WARNING:** The Frame should be well aligned with the Tomography's laser markings, to improve coordinate reliability.

10. MANUAL CALCULATION FOR STEREOTACTIC COORDINATES

Stereotactic calculations can be done manually, based on a computerized tomography exam.

WARNING: Find the exam image that best shows the target or region of interest.

Image example (Upward)

Figure 36

- **Determination of the “X” axis**

1. With the help of a Ruler, draw one line connecting points 1 and 3 and another connecting points 7 and 9;

Figure 37

2. Calculate and mark the mid-points (A and B) for each of these lines drawn;


Figure 38

3. Draw a straight line connecting points A and B; this line is, therefore, the X axis.


Figure 39

- **Determination of the “Y” axis**

1. With the help of a Ruler, calculate and mark the mid-point (C) of the X Axis (drawn in the previous step);


Figure 40

2. Draw a line connecting points 4 and 6;


Figure 41

3. Calculate and mark the mid-point (D) of the straight line drawn in the previous step;


Figure 42

4. Draw a straight line that passes through the mid-point (D) from the previous step and that passes through the mid-point (C) of the X axis. This line therefore represents the Y axis.


Figure 43

- **Determination of the “Z” coordinate**
- **For “UPWARD” target (positive Z Coordinate)**


Figure 44

- For targets in the RIGHT POSTERIOR quadrant, measure the distance between points 1 and 2 (Z_1);
- For targets in the LEFT POSTERIOR quadrant, measure the distance between points 8 and 9 (Z_3);
- For targets in the RIGHT and LEFT ANTERIOR quadrant, measure the distance between points 5 and 6 (Z_2);
- **For “DOWNWARD” target (negative Z Coordinate)**


Figure 45

- For targets in the RIGHT POSTERIOR quadrant, measure the distance between points 1 and 2 (Z_1);
- For targets in the LEFT POSTERIOR quadrant, measure the distance between points 8 and 9 (Z_3);
- For targets in the RIGHT and LEFT ANTERIOR quadrant, measure the distance between points 4 and 5 (Z_2);

- **Determination of the target's "X" and "Y" Coordinates**

The intersection of the lines drawn previously is the 0 point in "X" and 0 point in "Y." Therefore, based on these lines, we are able to locate the target coordinates, with the following steps:

1. With the help of a Ruler, draw a straight line parallel to the "X" axis, going through the target; draw a straight line parallel to the "Y" axis, going through the target


Figure 46

2. The intersection of these two new lines should be the target point (point E);


Figure 47

3. Measure the distance (X1) between point E and the "Y" axis. This distance is the "X" coordinate;


Figure 48

4. Measure the distance (Y1) between point (E) and the "X" axis. This distance is the "Y" coordinate;


Figure 49

NOTE: These coordinates correspond to the virtual coordinates, which should be converted into real coordinates using the method described in the next topic..

- **Conversion of virtual coordinates into real coordinates**

Print images are not always in real scale (1:1), since the properties can be changed to enter a greater number of images in a single print, causing reduction in these images; with this being the case, dimensions found in the print exam (which can be called virtual coordinates) should be processed and then entered into the equipment.

Virtual coordinates should be corrected using a conversion factor, which is found by dividing the Localization's Real Measurement (containing 140mm) by the Localization's Virtual Measurement, found by using a Ruler to measure the distance between points 1 and 3 in the print exam. As shown below:

$$\frac{140}{\text{MVF}} = \text{CONVERSION FACTOR}$$

After finding the conversion factor, the measurements for the virtual coordinates "X", "Y" and "Z" need to be multiplied by the value found, thus providing the real stereotactic coordinates.

- **Example:**

Suppose that the (virtual) measurements found in the CT exam were "X" = 22mm, "Y" = -18mm and "Z" = 18mm. And the measurement found between point 1 and point 3 was 73mm.

$$\frac{140}{73} = 1.92$$

The result (1.92) found is the conversion factor. Multiply the virtual coordinates found by the calculated conversion factor, according to the table below:

Axis	Virtual coordinate	Conversion factor	Stereotactic coordinate
X	22mm	x 1.92	42.24mm
Y	-18mm		-34.56mm
Z	18mm		34.56mm

The location of the target should be shown by the (+) and (-) minus signs:

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Z Coordinate: Top (+) and Bottom (-).

11. TRAYS FOR TRANSPORT / STORAGE / STERILIZATION

To transport and sterilize the Accessories, all parts must be placed in their appropriate locations.

-  **WARNING:** Store in a dry, ventilated and clean environment, protected from directed sunlight and substantial dust;
-  **WARNING:** Three boxes is the maximum amount recommended when piling for transport and storage;
-  **WARNING:** Before using the product, the package should be inspected as to its integrity, in order to guarantee that it is not compromised;
-  **WARNING:** Do not use the material if the package has been damaged and/or violated;
-  **WARNING:** Some of the information found on the package and product label include the maximum number of boxes that can be piled, protect from rain, fragile, and this side up.

12. CLEANING AND STERILIZATION

Cleaning method	
Sanitization	Wet a surgical compress or wipe with enzymatic detergent and apply to the entire surface of the item; rinse using a surgical compress or wipe moistened with water.
Manual Cleaning (immersion)	Immerse the dissembled parts in a container with a solution of cold to warm (up to 40°C) demineralized water and enzymatic detergent (dilute according to manufacturer recommendations). Let the solution act for at least 15 minutes. Using a soft brush, remove residue from the items, paying attention to recesses, screws and lights. Rinse in running demineralized water. Dry using a surgical compress or wipe and clean compressed air.
Automated Cleaning (Ultrasound)	Immerse the dissembled parts in an ultrasonic washing vat with a solution of cold to warm (up to 40°C) demineralized water and enzymatic detergent (dilute according to manufacturer recommendations). Activate the ultrasound for at least 15 minutes. Rinse in running demineralized water. Dry using a surgical compress or wipe and clean compressed air.
Sterilization Method	
SSUP / Autoclave	Sterilization in Saturated Steam under Pressure (SSUP / Autoclave): 134°C ±2 $\frac{2.14 \text{ kgf/cm}^2}{3 \text{ bar}}$ 7 minutes 100% H ₂ O dry for 50 minutes
ETO	Sterilization in Ethylene Oxide (ETO): 50°C ±10 $\frac{-0.3 \text{ kgf/cm}^2 \pm 0.1}{-0.3 \text{ bar} \pm 0.1}$ ≥40%UR 60 minutes 90%ETO and 10%CO ₂ 80 changes/hour ≥ 48 hours of aeration

The table below shows the cleaning and sterilization methods for the Stereotactic AimSystem TM03B Accessories:

Item	Code	Cleaning method			Sterilization Method	
		Sanitization	Manual Cleaning (Immersion)	Automatic Cleaning (Ultrasound)	Saturated Steam under Pressure (SSUP/ Autoclave)	Ethylene Oxide (ETO)
"A" Segment	TM0021FAA	✗	✓	✓	✓	✓
Ear Aligners	TM0004FAA	✗	✓	✓	✓	✓
Retractable Arms Support	OP1010FAA	✗	✓	✓	✓	✓
NMR Localizers	RM0002FAA	✓	✓	✓	✗	✗
Localizers for Magnetic Resonance	RM1000FAA	✓	✓	✓	✗	✗
Localizers for Computerized Tomography	CT0108FAA	✓	✗	✗	✗	✗
Support for Universal CT	CT100000FAA	✓	✗	✗	✗	✗
ELSCINT EXCEL 2400, 820	CT0800FAA	✓	✗	✗	✗	✗
GE, PACE, SYTE 1800, 3000, 400, SYNERGY PLUS	CT0200FAA	✓	✗	✗	✗	✗
GE9800	CT0300FAA	✓	✗	✗	✗	✗
GE HISPEED LXI, SX POWER, LIGHT SPEED ULTRA, LIGHT SPEED 8, 16, 32	CT1003FBA	✓	✗	✗	✗	✗
HITACHI W450, W4-40, W500, W2000, XR6000, PHILLIPS CX-S, CX-Q, SR4000, AV	CT1000FAA	✓	✗	✗	✗	✗
SHIMADZU	CT0400FAA	✓	✗	✗	✗	✗
PHILLIPS ALL MODELS UNSPECIFIED TOMOSCAN	CT0900FAA	✓	✗	✗	✗	✗
PHILLIPS CT16,40 BILLIANCE (UNTIL 11/2004)	CT-200000FAA	✓	✗	✗	✗	✗

PHILLIPS AURA	CT1004FAA	✓	✗	✗	✗	✗
TOSHIBA TCT300S, TCI300S, TCI500S, 600HQ, 600XY	CT0500FAA	✓	✗	✗	✗	✗
TOSHIBA AQUIL- ION/ASTEION	CT1005FAA	✓	✗	✗	✗	✗
PICKER PQ 2000, PQ 5000	CT1006FAA	✓	✗	✗	✗	✗
SIEMENS SOMAT- OM	CT0600FAA	✓	✗	✗	✗	✗
PHILLIPS BRILLI- ANE	CT201000FAA	✓	✗	✗	✗	✗
GE8800	CT0700FAA	✓	✗	✗	✗	✗

Key			
✓	Procedure Permitted	✗	Procedure Prohibited

⚠ WARNING: Items exposed to Creutzfeldt-Jakob disease should be discarded according to local regulations.

⚠ WARNING: If the parameters for cleaning and sterilization are not followed as described here, contamination control may not be guaranteed or damages could be caused to the product.

⚠ WARNING: Do not use abrasive and corrosive agents for cleaning, since they can cause damage to the product.

⚠ WARNING: After cleaning, the components should be immediately dried using a dry compress and compressed air. Pay attention to recesses and screws.

NOTE: Follow the washing machine manufacturer's instructions.

NOTE: Follow the information specified by the detergent manufacturer.

13. LUBRICATION

Prior to the surgery, use sterilized vaseline or silicone to lubricate sliding parts and the tips of fixation screws. Lubrication of each part makes tightening and loosening screws easier and facilitates movement of parts during surgery.

14. PREVENTIVE AND CORRECTIVE MAINTENANCE

To guarantee the appropriate functioning and performance of the equipment, it should be returned for preventive maintenance after 18 months.

If there is any defect or malfunction before 18 months have passed, do not use the equipment and immediately contact Micromar or your authorized sales representative.

Micromar is solely authorized to provide technical assistance for the Stereotactic AimSystem TM03B. After the 18 month period has passed, the equipment must be sent every 14 months to Micromar for assessment to find whether any maintenance or disposal is needed.

PREVENTIVE MAINTENANCE TABLE			
Quality Control	Quality Control	Quality Control	Quality Control
1 st Preventive Maintenance	2 nd Preventive Maintenance	3 rd Preventive Maintenance	4 th Preventive Maintenance

15. SYMBOLS

Symbol	Key
	Not sterile
	Keep out of sunlight
	Manufacturer
	Attention, check documents
	Keep out of rain
	Operating manual, operating instructions
	Date of manufacture
	Serial number
	Stacking limit by number
	Fragile
	This side up
	Lot number

16. CARE AND PRECAUTIONS

- The equipment should only be used with the basic component items containing the same serial number (except accessories). Replacement of components with components containing other serial numbers could compromise the equipment's precision, causing damage to the patient.
- Accidental movements involving the accessory components of the Stereotactic AimSystem TM03B during the procedure can compromise the efficacy of the procedure.
- Avoid dropping, twisting, crushing, impacting or any action that will deform accessories.
- In the event of any of the above occurrences, the user should contact Micromar for assessment of the equipment.
- Deformation of the accessory structure compromises the system's precision, placing the patient's life at risk.
- Incorrect reading of accessory scales (user error in reading or because of poor equipment condition) can put the patient's life at risk.
- Screw in and correctly attach the movable parts of the accessory.
- To perform the procedure, the accessory should be cleaned and sterilized according to the instructions in this manual.
- Immediately after the end of the procedure, the equipment should be cleaned according to the instructions in this manual.
- The accessory should be sterilized according to the instructions contained in this User Manual; otherwise, there is a risk of bacterial or viral contamination, causing infectious complications for the patient.
- This product may only be manipulated by a surgeon specialized in neurosurgery and trained by Micromar.
- The surgeon must consider the patient's physical and pathological conditions when attaching the Stereotactic AimSystem TM03B fixation screw to the patient's head, in order to prevent damage to same.
- If the instructions contained in the manual and in the user training are not followed, this could result in incorrect use of the System, placing the patient's life at risk.
- Use the maintenance service provided by Micromar or an authorized representative to correctly repair the system if any deformations occur.
- To repair or replace any component of this accessory, you must contact Micromar through customer service or through any authorized representative.
- Any accessory used to treat patients suffering from Creutzfeldt-Jakob disease MUST NOT BE reused for other patients. Check local disposal laws.
- The Stereotactic AimSystem TM03B is solely meant for use in hospitals and clinics.
- Adverse effects: There are no secondary or side effects directly related to appropriate use of the product.

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QUALITY CERTIFICATION

This manual is aimed at helping the user understand how to correctly, safely and effectively use the Stereotactic Aymsystem TM03B. The user must carefully read the user manual that comes with the product before using it.

This product is warranted by Micromar Ind. e Com. Ltda. within the following conditions:

- 1- An 18 month warranty, as from the date of issue of the receipt of sale to the consumer, for replacement of parts or components that may have manufacturing defects, except for those resulting from improper use or installation, done not following the specifications contained in the user manual;
- 2- Warranty pending presentation of this certificate will be offered at no cost to the customer;
- 3- This warranty is null and void in the event of:
 - 3.1- Defect caused by the customer or by third parties other than the manufacturer;
 - 3.2- The product has been changed or repaired by persons or companies not authorized by the manufacturer;
 - 3.3- The recommendations of the user manual are not followed;

Serial Number: _____

Invoice Number: _____

Date: _____

EQUIPMENT INSPECTED BY:

Quality Control

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
"A" Segment	TM002IFAA	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
Localizers for nuclear magnetic resonance	RM1000FAA RM0002FAA	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
Ear Aligners	TM0004FAA	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
Retractable Arms Support	OP1010FAA	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
Isocentric retractor spatulas	OP0010FAA	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

Product	Code	SN	Date	PREVENT MAINTENANCE TABLE			
Support for CT _____	_____	SN: _____	__/__/__	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL	__/__/__ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

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