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
M I C R O M A R

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FIXADOR DE CRÂNIO HEADFIX HF03B

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PRODUCT DESCRIPTION:

HEADFIX HF03B CRANIAL FIXATION

ONLY USE THESE ACCESSORIES AFTER READING THE INSTRUCTIONS AND RESPECTIVE WARNINGS

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MICROMAR IND. COM. RESERVES THE RIGHT TO MAKE CHANGES TO THIS MANUAL WITHOUT PRIOR NOTICE.

CONTACT OUR SALES DEPARTMENT TO RECEIVE THE LATEST VERSION OF THE USER MANUAL.

THIS MANUAL CONTAINS INFORMATION PROTECTED BY COPYRIGHT.

THIS MANUAL MAY NOT BE REPRODUCED, TRANSLATED OR TRANSMITTED BY ANY ELECTRONIC OR OTHER MEANS, IN WHOLE OR IN PART, WITHOUT EXPRESS AUTHORIZATION FROM MICROMAR.

INDEX

1. INTRODUCTION.....	30
2. ILLUSTRATIONS.....	30
3. USE OF ACCESSORIES.....	30
4. FRAME ADAPTOR - HF5000FAA	31
5. DOUBLE FIXER - HF6000FAA	33
6. SIDE CLAMP - HF7000FAA.....	35
7. HEADREST (HORSESHOE) - HF8000FAA	38
8. SET OF 2 GEL PADS - LEFT AND RIGHT - MF0038FAA.....	41
9. CROSSBAR ADAPTOR FOR SEATED POSITION - HF11000FAA.....	43
10. TRANSPORT AND STORAGE CASE - BASIC.....	46
11. TRANSPORT AND STORAGE CASE - BASIC AND OPTIONAL	49
12. CLEANING AND STERILIZATION.....	54
13. PREVENTIVE AND CORRECTIVE MAINTENANCE	56
14. SYMBOLS.....	57
15. TRANSPORT AND STORAGE	57
16. WARNINGS AND PRECAUTIONS.....	58

17. INTRODUCTION

This manual is aimed at helping the user of the HeadFix HF03B Cranial Fixation System and its accessories to safely and effectively use the equipment.

Before using the accessories, the user should carefully read the manual and gain a total understanding of its contents.

This manual should be kept with the equipment for reference in the event that questions arise.

18. ILLUSTRATIONS

The images in this manual show an illustration of the accessories. There may be some differences in comparison to the physical components. These differences in the illustrations do not interfere with the use of the accessories.

19. USE OF ACCESSORIES

The accessories must be assembled and used as indicated in this manual and for the purpose for which they were created.

The accessories described herein can only be used in combination with other equipment or components that are expressly recognized by Manufacturer as being compatible therewith.

Compatible products:

- Heafix HF03B cranial fixation system;
- Micromar Stereotactic TM03B;
- HeadFix HF02B Cranial Fixation;
- AimNav instruments - Cranial Toolbox - NN00100FAA.

20. FRAME ADAPTOR - HF5000FAA

The Frame adapter allows for connection of the Cranial Fixation with the Micromar Frame.

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.


Figure 1

1. Turn the handle counterclockwise (⚙️) **I** and then leave it in a position that allows for insertion of the frame **II**;


Figure 2

2. Insert the threaded pin for the Articulated system into the connection on the Frame Adaptor and adjust its position until the teeth fit together uniformly (I), then turn the Articulated system handle clockwise (II) to attach the Adapter to the equipment (II);


Figure 3

3. Insert the Frame in the Frame Adaptor connection profile;


Figure 4

4. Position the lock over the Frame (I), then turn the handle clockwise (II) to attach the Frame to the Frame Adaptor (II);


Figure 5

⚠ WARNING: Ensure the correct placement of the fixation teeth during fixation;

⚠ WARNING: See chapter 12 for more information on the correct cleaning process;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

21. DOUBLE FIXER - HF6000FAA

Double fixer is an accessory that allows for use of Micromar navigation articulated arms in conjunction with the cranial fixation system during the surgical procedure.


Figure 6

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.

1. Insert the threaded pin for the Double fixer into the connection of the Skull clamp and adjust its position until the teeth fit together uniformly (I), then turn the handle on the Double fixer clockwise (II) to attach it to the equipment (II):


Figure 7

2. Insert the Micromar navigation articulated arm threaded pin for into the Double fixer connection and adjust its position so that the teeth fit together uniformly. Next, turn the handle on the articulated arm clockwise to attach;


Figure 8

⚠ WARNING: Ensure the correct placement of the fixation teeth during fixation;

⚠ WARNING: See chapter 12 for more information on the correct cleaning process;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

22. SIDE CLAMP - HF7000FAA

The side clamp is an accessory that allows the HF03B Headfix Cranial Fixation System to be used with Leyla surgical retractors.


Figure 9

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.

3. Send the side clamp with your sterilization and storage case to the material sterilization center;
4. After the sterilization process, remove the cover from the case pressing the side latches **I**, then gently pull the cover **II**;


Figure 10

5. Remove the side clamp from the case;


Figure 11

6. Turn the handle counterclockwise (🔓) until it's possible to couple the side clamp on the Skull Clamp;


Figure 12

7. Insert the side clamp over the sterile field on the Skull Clamp (I), then turn the handle clockwise to (🔒) para to attach the side clamp to the skull clamp (II);


Figure 13

8. Turn the second handle counterclockwise (🔓) until it's possible to couple the Leyla retractor arm to the side clamp fitting (I), then couple the retractor arm to the side clamp (II);


Figure 14

9. Turn the handle clockwise () to attach the retractor arm to the side clamp.


Figure 15

 **WARNING:** Ensure the correct attachment of the side clamp to the HF03B Cranial Fixation System during the assembly;

 **WARNING:** Ensure the correct attachment of the retractor arm to the side clamp during the assembly;

 **WARNING:** Ensure the correct placement of the fixation teeth during fixation;

 **WARNING:** See chapter 12 for more information on the correct cleaning process;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

23. HEADREST (HORSESHOE) - HF8000FAA

The Headrest (Horseshoe) is recommended for procedures where it's necessary to support and stabilize the patient's head during the surgical procedure in the prone position.


Figure 16

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment. [RM#MANUAL#Q16]

1. Insert the threaded pin for the Articulated system into the connection on the Headrest and adjust its position until the teeth fit together uniformly **I**, then turn the handle clockwise **II** to attach the Headrest to the equipment;


Figure 17

2. With all of the equipment's articulations unlocked, adjust the position according to the need of the surgical procedure, then lock all articulations;


Figure 18

3. Insert the set of 2 Gel Pads;


Figure 19

4. Turn the handle located on the bottom of the Headrest counterclockwise (🔓) to release the adjustment movement for the headrests (I) and then adjust the spaces between them according to the needs of the procedure (II);


Figure 20

5. After adjustment, turn the handle clockwise (🔒) to lock the adjustment movement of the headrests (I), then position the patient's head on the set of 2 Gel Pads (II);


Figure 21

• **Preparation for assembly of the cervical traction rod**

6. Turn the handle on the cervical traction rod counterclockwise (🔓) to release the adjustment movement (I) after adjusting the rod's position according to the needs of the procedure (II);


Figure 22

7. Turn the handle clockwise (🔒) to attach the rod to the Headrest;


Figure 23

⚠ WARNING: Ensure the correct placement of the fixation teeth during fixation;

⚠ WARNING: See chapter 12 for more information on the correct cleaning process;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

NOTE: See Chapter 8 for more information on how to use the set of Gel Pads.

24. SET OF 2 GEL PADS - LEFT AND RIGHT - MF0038FAA

The Set of 2 Gel Pads is recommended for use in conjunction with the Headrest (Horseshoe), to support and stabilize the patient's head during the surgical procedure in the prone position.


Figure 24

1. Check that the Set of 2 Gel Pads contains fastening clips;


Figure 25

2. Position the fastening clips over the left and right support holes on the headrest **I**, then press them until the sets connect **II**;


Figure 26

3. After the end of the procedure, hold the plastic part of the set of 2 Gel Pads **I** and then pull the set to disassemble it **II**.


Figure 27

⚠ WARNING: See chapter 12 for more information on the correct cleaning process;

⚠ WARNING: The set of 2 Gel Pads must not be sterilized;

⚠ WARNING: Sharp objects may damage the set of 2 Gel Pads;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

NOTE: See Chapter 7 for more information on how to use the set of Headrest (Horseshoe).

25. CROSSBAR ADAPTOR FOR SEATED POSITION - HF11000FAA

The Crossbar adaptor for seated position helps the surgeon to perform procedures in the sitting position


Figure 28

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.

1. Remove the handle from the tube, turning it counterclockwise **I** and then remove it the Base **II**;


Figure 29

2. After removing the Base's handle, first lift the hex key **I** until it's perpendicular to the handle **II**;


Figure 30

3. Insert the hex key into the crossbar adapter screws **I**, then turn the handle counterclockwise until the rods can be accessed **II**;


Figure 31

4. Insert the rods into the holes **I**, with the hex key inserted into the screw, turn the handle clockwise to attach them **II**;


Figure 32

5. With the patient positioned, insert the crossbar adapter for seated position rods into the fittings on the surgical table;


Figure 33

6. Lock the surgical table connection adapters to attach the crossbar adapter for seated position;
7. Insert the hex key into the Crossbar adaptor's screws (I), then turn the handle counterclockwise to release access to the rods on the Base (II);


Figure 34

8. Insert the Base set's rods into the holes for the crossbar adapter for seated position (I), then turn the cross-bar adapter screws clockwise to attach the equipment (II);


Figure 35

⚠ WARNING: See chapter 12 for more information on the correct cleaning process;

NOTE: For more information on the correct use of the HF03B Cranial Fixation System, consult the user manual for said equipment.

26. TRANSPORT AND STORAGE CASE - BASIC

The case assists in the transportation and storage of the equipment.


Figure 36

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.

1. Pull the side latches on the case to unlock the cover **I**, then open the case **II**;


Figure 37

2. Remove the mobile support using the strap;


Figure 38

3. Remove the base set;


Figure 39

4. Remove the articulated system set;


Figure 40

5. Remove the Skull Clamp set;


Figure 41

6. After using and cleaning of the equipment, place the sets their proper places;


Figure 42

7. Insert the mobile support;


Figure 43

8. Close the cover **I**, then lock the case **II**;


Figure 44

⚠ WARNING: Make sure the case is locked before transport it.

27. TRANSPORT AND STORAGE CASE - BASIC AND OPTIONAL

The case assists in the transportation and storage of the equipment and its accessories.


Figure 45

⚠ WARNING: For more information on the associated equipment used, consult the user manual for said equipment.

1. Pull the side latches on the case to unlock the cover **I**, then open the case **II**;


Figure 46

2. Remove the base set;


Figure 47

3. Remove the articulated system set;


Figure 48

4. Remove the skull clamp set;


Figure 49

5. To remove the accessories from case it's necessary to remove the mobile support using the straps;


Figure 50

6. Remove the frame adaptor;


Figure 51

7. Remove the double fixer;


Figure 52

8. Headrest (Horseshoe) with the gel supports;


Figure 53

9. Remove the crossbar adaptor with its rods.


Figure 54

10. After using and cleaning of the accessories, place the sets their proper places;


Figure 55

11. Insert the mobile support;


Figure 56

12. Place the equipment sets in their proper places;


Figure 57

13. Close the cover **I**, then lock the case **II**;


Figure 58

⚠ WARNING: Make sure the case is locked before transport it.

28. CLEANING AND STERILIZATION

The table below shows the methods of cleaning and sterilizing the accessories:

Cleaning Method	
Sanitation	Wet a surgical compress or wiper with enzymatic detergent and apply over the entire surface of the item; rinse with surgical compress or wiper moistened with water.
Manual Cleaning (immersion)	Immerse the disassembled parts in a container with cold to warm demineralized water (up to 40°C) with enzymatic detergent (dilution as recommended by the manufacturer). Allow the solution to act for at least 15 minutes. With a soft brush remove the residue from the items, taking care of the recesses, threads and lumens. Rinse with running demineralized water. Dry with surgical compress or wipe and clean compressed air.
Cleaning Automated (Ultrasound)	Immerse the disassembled parts in an ultrasonic washer tank with cold to warm demineralized water solution (up to 40°C) with enzymatic detergent (dilution as recommended by the manufacturer). Operate the ultrasound for at least 15 minutes. Rinse with running demineralized water. Dry with surgical compress or wiper and clean compressed air.
Sterilization Method	
Steam	Sterilization in Saturated Vapor under Pressure (Steam): 134°C ±2 $\frac{2.14 \text{ kgf/cm}^2}{3 \text{ bar}}$ 7 minutes 100% H ₂ O 50 minutes of drying
ETO	Sterilization in Ethylene Oxide (ETO): 50°C ±10 $\frac{-0.3 \text{ bar} \pm 0.1}{-0.3 \text{ kgf/cm}^2 \pm 0.1}$ ≥40%UR 60 minutes 90%ETO e 10%CO ₂ 80 exchanges/hour ≥ 48 hours of aeration

Item	Code	Cleaning Method			Sterilization Method	
		Sanitation	Manual Cleaning (immersion)	Cleaning Automated (Ultrasound)	Saturated Steam Under Pressure (Steam)	Ethylene Oxide (ETO)
Frame adaptor	HF5000FAA	✗	✓	✓	✗	✗
Double fixer	HF6000FAA	✗	✓	✓		
Headrest (Horseshoe)	HF8000FAA	✗	✓	✓		
Set of 2 Gel Pads - left and right	MF0038FAA	✓	✗	✗		
Crossbar Adaptor for Seated Position	HF11000FAA	✗	✓	✓		
Side clamp	HF7000FAA	✗	✓	✓	✓	✓
Subtitle		✓	Allowed Procedure		✗	Prohibited Procedure

-  **WARNING:** Items exposed to Creutzfeldt-Jakob disease should be discarded according to local regulations;
-  **WARNING:** If the parameters for cleaning and sterilization are not followed as described herein the contamination control may not be guaranteed or damage may be caused to the product;
-  **WARNING:** Do not use abrasive and corrosive agents for cleaning, since they can cause damage to the product;
-  **WARNING:** After cleaning, the components should be immediately dried using a dry compress and compressed air. Pay attention to recesses and screws;

NOTE: Follow the manufacturer instructions for the washing machine;

NOTE: Follow the information specified by the detergent manufacturer.

29. PREVENTIVE AND CORRECTIVE MAINTENANCE

To guarantee the appropriate functioning and performance of the accessories, they must return for preventive maintenance **after 18 months**.

If there is any defect or malfunction before 18 months have passed, do not use the accessory and immediately contact Micromar or your authorized sales representative.

Micromar is the only authorized technical service provider for the Headfix HF03B Cranial Fixation System accessory line.

After the expiration of the 18-month period, it's need sent the accessories to Micromar, and the **every 14 months**, for assessment and to perform any maintenance or disposal.

PREVENTIVE MAINTENANCE TABLE			
____/____/____ Quality control	____/____/____ Quality control	____/____/____ Quality control	____/____/____ Quality control
1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

30. SYMBOLS

Symbols	Subtitle
	Keep away from sunlight/Keep away from heat
	Manufacturer
	Attention: Consult the instructions for use
	Keep dry/Keep away from humidity
	Authorized representative in the European Community
	Consult the instructions for use
	Date of manufacture
	Serial number
	Stacking limit
	Fragile
	This way up
	Lot number

31. TRANSPORT AND STORAGE

The product must be transported and stored in the packaging provided by Micromar.

32. WARNINGS AND PRECAUTIONS

- Accidental movements involving accessories during the procedure can compromise the efficacy of the procedure.
- Avoid dropping, twisting, crushing, impacting or any action that will deform the accessories.
- In the event that one of the above occurs, the user should contact Micromar for assessment.
- Screw in and correctly attach the movable parts of the accessories.
- To perform the procedure, the accessories should be cleaned according to the instructions in this manual
- Immediately after ending the procedure, the accessories must be cleaned according to the instructions in this manual.
- These accessories are solely made for use in surgical centers.
- These accessories may only be manipulated by a surgeon specialized in neurosurgery with experience performing cranial fixation.
- If the instructions contained in the manual user supplied by Micromar are not followed, this could result in incorrect use of the accessories, placing the patient's life at risk.
- Use the maintenance service provided by Micromar or an authorized representative to correctly repair the accessories when there are any problems with its functioning.
- To repair or replace any component of these accessories, you must contact Micromar through customer service or through any authorized representative.
- Accessories used to treat patients suffering from Creutzfeldt-Jakob disease MUST NOT BE reused for other patients.
- Disposal of the product should be done according to local laws.
- The accessories have a useful life of 5 years, as from the date of issuance of the invoice.

Avenida Marginal ao Córrego da Serraria, nº168 • SP • Brazil • Phone: (+55 11) 4057-5722

QUALITY CERTIFICATION

CRANIAL FIXATION SYSTEM AND ACCESSORIES

This product is guaranteed by Micromar Ind. e Com. Ltda. under the following conditions:

- 1- An 18 month warranty, as from the date of issue of the receipt of sale to the consumer, for replacement of parts or components that may have manufacturing defects, except for those resulting from improper use or installation, done not following the specifications contained in the user manual;
- 2- Warranty pending presentation of this certificate will be offered at no cost to the customer;
- 3- This warranty will be null in the event:
 - 3.1- Of defect caused by the customer or by third parties other than the manufacturer;
 - 3.2- The product has been changed or repaired by persons or companies not authorized by the manufacturer;
 - 3.3- The recommendations of the user manual are not followed;
- 4- EAny defects resulting from natural wear to the product or caused by agents of nature are not covered by this warranty. After expiration of the stipulated term, the HeadFix and Accessories should be sent to Micromar, for assessment of any maintenance and/or disposal of the product.
- 5- After the warranty has expiry and the every 14 months, the Cranial fixation system and the accessories must be sent to Micromar for evaluation of the equipment to detect possible maintenance or disposal of the equipment.

*** This manual and certificate should be kept with the product and they may be sent along with the equipment each time it is sent for preventive or corrective maintenance.**

Equipment inspected by:

Quality Control

Product	Code	SN	Date	PREVENTIVE MAINTENANCE TABLE			
Frame adaptor	HF5000FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance
Double fixer	HF6000FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance
Side clamp	HF7000FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance
Headrest (Horseshoe)	HF8000FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance
Set of 2 Gel Pads - left and right	MF0038FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance
Crossbar Adaptor for Seated Position	HF11000FAA	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

Product	Code	SN	Date	PREVENTIVE MAINTENANCE TABLE			
Transport and Storage Case - Basic	HF15000FAA	SN.:_____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

Product	Code	SN	Date	PREVENTIVE MAINTENANCE TABLE			
Transport and Storage Case - Basic and Optional	HF16000FAA	SN.:_____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

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