DoubleArc®

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Português · English · Español

MICROMAR

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ESTEREOTÁXICO DOUBLEARC TM04B

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PRODUCT DESCRIPTION:

STEREOTACTIC DOUBLEARC TM04B

ONLY USE THE EQUIPMENT AFTER CAREFULLY READING THE USER'S MANUAL AND RESPECTIVE WARNINGS.
CHECK IF THIS VERSION OF INSTRUCTION FOR USE IS ACCORDING TO THE PRODUCT PURCHASED.

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INFO@MICROMAR.COM

WWW.MICROMAR.COM

MICROMAR IND. COM. THE RIGHT TO MAKE CHANGES TO THIS MANUAL WITHOUT PRIOR NOTICE IS HEREBY RESERVED. CONTACT OUR SALES DEPARTMENT TO RECEIVE THE LATEST VERSION OF THE USER MANUAL.

THIS MANUAL CONTAINS INFORMATION PROTECTED BY COPYRIGHT

THIS MANUAL MAY NOT BE REPRODUCED, TRANSLATED OR TRANSMITTED BY ANY ELECTRONIC OR OTHER MEANS, IN WHOLE OR IN PART, WITHOUT EXPRESS AUTHORIZATION FROM MICROMAR

CONTENTS

1. INTRODUCTION.....	46
2. ILLUSTRATIONS.....	46
3. COMPATIBILITY.....	46
4. USE OF EQUIPMENT	46
5. PRESENTATION	46
6. BASIC CONCEPTS AND OPERATION PRINCIPLE	47
7. DOUBLEARC TM04B	48
8. INDICATIONS FOR USE	48
9. PROCEDURES	48
9.1. COMPONENT LIST.....	49
9.2. IDENTIFICATION OF THE TRAYS	52
9.3. PREPARING THE PATIENT FOR THE SURGICAL PROCEDURE - FIXATION TRAY.....	53
9.3.1. DIRECT FIXATION.....	54
9.3.2. INDIRECT FIXATION.....	56
9.3.3. PREPARATION FOR IMAGING EXAM	59
9.3.4. MANUAL CALCULATION OF STEREOTACTIC COORDINATES.....	60
9.3.5. ASSEMBLY OF THE COMPONENTS - COORDINATE TRAYS.....	68
10. SCALES.....	77
11. TRAYS FOR TRANSPORT / STORAGE / STERILIZATION.....	77
12. CLEANING AND STERILIZATION.....	78
13. LUBRICATION.....	79
14. PREVENTIVE AND CORRECTIVE MAINTENANCE.....	80
15. SYMBOLS	81
16. WARNINGS AND PRECAUTIONS.....	82

1. INTRODUCTION

This manual is intended to instruct the user on the safe and effective use of the DoubleArc TM04B Stereotactic System.

Before using the product, the user must read the manual thoroughly, fully understand its contents, and complete the training provided by Micromar.

This manual should be kept with the equipment for reference in case of any questions.

2. ILLUSTRATIONS

The images in this manual illustrate the product and its accessories. There may be differences compared to prior equipment models and versions. These differences in illustrations do not interfere in its use.

3. COMPATIBILITY

The equipment described in this manual may only be used in combination with other components that are expressly recognized by Micromar as compatible.

Compatible products:

- Micromar TM03B Stereotactic System accessories;
- Micromar Headfix Frame Adapter;
- Micromar Stereotaxy Kits;

For further information, please contact Micromar.

4. USE OF EQUIPMENT

The DoubleArc TM04B Stereotactic System is intended to be assembled and used only in accordance with the procedures described in this manual and for the purpose for which it was developed.

5. PRESENTATION

The DoubleArc TM04B Stereotactic System is supplied in its basic form, containing all components necessary to perform the surgical procedures for which the product was developed.

6. BASIC CONCEPTS AND OPERATION PRINCIPLE

Etymologically, stereotaxy (from Greek, stereos - solid, three-dimensional; taxis - placement, position) means guided three-dimensional localization. Stereotactic neurosurgery means reaching a target on a region of the brain, using Cartesian coordinates that are defined in advance using a geometric orientation device, which is attached to the head and allows for an instrument to be guided to a target.

The Stereotactic DoubleArc TM04B allows for targets to be located using Cartesian coordinates (X, Y and Z), as shown in Figure 1.


Figure 1

"X" = Right and left sides

"Y" = Back and front

"Z" = Top and bottom

The equipment uses an isocentric concept that allows a target to be reached from any point of entry; in other words, the target is reached with precision regardless of Alpha and Beta angles. Figure 2 shows how the equipment can be angled.


Figure 2

7. DOUBLEARC TM04B

Technical information of the equipment:

Net weight:

- 3,5kg;

Height:

- 400mm;

Width:

- 593 mm;

Isocenter:

- 243 mm;

Precision:

- Linear precision: 0.5 mm;
- Angular precision: 1°.

8. INDICATIONS FOR USE

The Micromar DoubleArc TM04B Stereotactic System enables intracranial access through neurosurgical procedures.

1. In preparing the procedure, stabilizing the patient's position for imaging exams, such as Magnetic Resonance and Computerized Tomography..
2. Targets' locations obtained through the patient's imaging exam.

9. PROCEDURES

The equipment was developed to support and assist in the performance of the following procedures:

Guided Neurosurgery:

- Guided surgery;
- Catheter implantation;
- Guided neuroendoscopy;;

Diagnostic Neurosurgery:

- Biopsy;

Therapeutic Neurosurgery:

- Haematoma Evacuation;
- Abscess Evacuation;

Functional Neurosurgery:

- Treatment of movement disorders (Lesioning and Deep electrode-DBS Implantation);
- Treatment of chronic pain
- Behavioral Surgery
- Spasticity
- Treatment of Epilepsy

9.1. COMPONENT LIST

The figures below show the composition of the DoubleArc TM04B Stereotactic System model.

• **Fixation Tray**

The fixation tray is composed of:

Item	Name	Code	Quantity
1	Sterilization Tray 1	TM0030CAA	1
2	Frame	TM0020FAA	1
3	Reusable Fixation Screws (Stainless Steel)	TM0230FAA	1
4	Extension Posts	TM0140FAA	1
5	Allen Key 4.0mm	TM0032FAA	2
6	Allen Key 6.0mm	TM0034FAA	2
7	Curved Segment	TM0020FAD	1


Figure 3

• **Reusable Fixation Screws (White Polyacetal)**

Item	Name	Code	Quantity
8	Reusable Fixation Screws (White Polyacetal)	TM0126FAA	1


Figure 4

• **Coordinate Tray A**

The Coordinate Tray A is composed of:

Item	Name	Code	Quantity
9	Sterilization Tray 2	TM0040CAA	1
10	Millimetrical Ruler*	TM0024FAA	1
11	Fixation Screw	PA0020FAA	3
12	“Z” Guide Post*	TM0011FAA	1
13	“Z” Post*	TM0113FAA	1
14	Short “Z” Guide Post*	TM0012FAA	1
15	Short “Z” Post*	TM0114FAA	1
16	Millimetrical Tube*	TM0022FAA	1
17	Arc*	TM0025FAA	1
18	Chariot*	TM0010FAA	1

• **Coordinate Tray B**

The Coordinate Tray B is composed of:

Item	Name	Code	Quantity
19	Sterilization Tray 2	TM0040CAA	1
20	Millimetrical Ruler*	TM0024FAA	1
21	Fixation Screw	PA0020FAA	3
22	“Z” Guide Post*	TM0011FAA	1
23	“Z” Post*	TM0113FAA	1
24	Short “Z” Guide Post*	TM0012FAA	1
25	Short “Z” Post*	TM0114FAA	1
26	Millimetrical Tube*	TM0022FAA	1
27	Arc*	TM0025FAA	1
28	Chariot*	TM0010FAA	1


Figure 5

NOTE: The items marked with * on Coordinate Trays A and B are identified by the corresponding letters A and B.

Item	Name	Code	Quantity
29	Sight	TM1022FAA	1


Figure 6

NOTE: The DoubleArc TM04B Stereotactic System has only one sight, which may be placed either in the graduated tube A or in the graduated tube B.

• **Support for Basic Computerized Tomography**

Item	Name	Code	Quantity
19	Support for Basic CT	CT0100FAA	1


Figure 7

9.2. IDENTIFICATION OF THE TRAYS

The sterilization trays have specific features that allow for their correct identification even after sterilization and when covered by the surgical drape.

- **Fixation Tray**

The Fixation Tray is identified by the presence of the frame and the absence of support feet, as illustrated in the figure below.


Figure 8

- **Coordinate Trays A and B**

The Coordinate Trays A and B are identified by the presence of the arch and the four support feet, making it possible to recognize them by touch when sterilized and covered by the surgical drape, as illustrated in the figure below.


Figure 9

NOTE: For information on proper sterilization, refer to Chapter 12.

9.3. PREPARING THE PATIENT FOR THE SURGICAL PROCEDURE - FIXATION TRAY

To correctly attach the Frame to the patient, the position of the target needs to be verified using a previous imaging exam study.

The Frame should be positioned above or below the target, never in the region of interest, since this will prevent it from being viewed on magnetic resonance and/or computerized tomography exams.

Attachment of the Frame to the patient's head can basically occur in two ways, using DIRECT FIXATION, with the Frame attached to the skull by screws that are screwed directly to the Frame, or INDIRECT FIXATION, using Extension Posts (Optional) that attach the Frame to the patients head using screws..

- **Place the Frame on the patient's head**

Figures (a) and (b) show how accurately the Frame is positioned to perform the surgical procedure.

When the target to be used is located above the Frame, the UPWARD (a) procedure should be followed; and when the target to be used is below the Frame, the DOWNWARD (b) procedure should be followed.


Figure 10

⚠ WARNING: Micromar recommends that the Frame be positioned ± 2 cm above or below the target to be used during the surgical procedure;

⚠ WARNING: The Frame should never be positioned at the same height as the region of interest;

9.3.1. DIRECT FIXATION

DIRECT FIXATION of the Frame on the patient's head can occur using 3 or 4 fixation points, as follows:

- **Direct fixation of the Frame using three screws**

For direct fixation using three screws, the screws are screwed directly into the Frame using the Straight Segment.

If the surgical approach is made using the Ruler attached to the front region of the Frame, a medium or small screw needs to be used on the front part (depending on the size of the patient's head), using two medium or large screws on the back part. Care should be taken so that the front screw does not extend beyond 3 mm from the face (A) of the Frame; otherwise, it will prevent the Ruler from being appropriately placed on the front region.


Figure 11

⚠ WARNING: Fixation screws should not extend beyond faces (A) and (B);

⚠ WARNING: Fixation screws have a piercing tip and should therefore be handled with care.

- **Direct fixation of the Frame using four screws**

For direct fixation using four screws, two front screws and two back screws will be used. Screws should not extend beyond faces (A) and (B) of the Frame, since this prevents the Ruler from being attached correctly.


Figure 12

⚠ WARNING: Fixation screws should not surpass faces (A) and (B);

⚠ WARNING: Fixation screws have a piercing tip and should therefore be handled with care.

9.3.2. INDIRECT FIXATION

The Stereotactic DoubleArc TM04B also allows for indirect fixation that allows the Frame to be positioned lower on the skull, which allows targets located at the base of the brain to be reached.

- **Indirect fixation of the Frame using four screws**

For indirect fixation using four screws, the Curved Segment and Extension Posts will be used (in addition to the Frame):

- CURVED SEGMENT

Recommended use

The Curved Segment is recommended when a patient needs to be intubated and the Frame is directly or indirectly attached with four screws. The need to intubate the patient should be assessed in advance, since the system's assembly can not be changed during the surgical procedure.


Figure 13

Curved Segment Assembly:

1. To assemble the Curved Segment on the Frame, the Straight Segment must first be removed;
2. After removing the Straight Segment, insert the Curved Segment so that it easily fits;
3. After placing the Curved Segment, tighten the screws (the same ones removed from the Straight Segment), turning them clockwise by using a 4mm Allen Key;
4. To remove the Curved Segment, release the screws holding the segment by using a 4mm Allen Key, turning it counterclockwise.


Figure 14

⚠ WARNING: Never change the segment when the Frame is attached to the patient's head, since this will cause irreversible damage the Frame;;

⚠ WARNING: If the Curved Segment needs to be positioned downward, turn the Frame assembled to the Curved Segment before attaching it to the patient's head;

- EXTENSION POSTS

Recommended use

The Extension Posts are used for surgical approaches (always Upward), with indirect fixation, where the target is found on the more basal region of the brain; they are also used in procedures where a greater brain region range is desired.


Figure 15

Assembly of the Extension Posts

1. Remove the screws from the Extension Posts, turning them counterclockwise;


Figure 16

2. Insert the Extension Posts, each in their position, in the Frame's internal part (Figure 16);
3. Lock the Frame posts using the screws (removed in the first step) and use the 4mm Allen Key to attach the Extension Posts to the Frame;


Figure 17

4. If the Reusable Stainless Steel or Disposable Aluminum Screws are used, then the Adapter Bushing Guides need to be used (in this case, use the 4mm Allen Key to tighten the screws);


Figure 18

5. To assemble the Adapter Bushing Guides, they must be turned clockwise on the inside of the Frame and the screws should be inserted from the outside to the inside.


Figure 19

6. If the Reusable Polyacetal Screws are used, remove the threaded Bushing Guides need to be used (in this case, use the 6mm Allen Key to tighten the screws).

Figure 20

WARNING: Fixation Screws have a piercing tip and should therefore be handled with care.

9.3.3. PREPARATION FOR IMAGING EXAM

Micromar has accessories to help the user in obtaining images on equipment for:

- Localizers for Computerized Tomography


Figure 21

- The set consists of a frontal reference, used in the frontal quadrant of the stereotactic frame, and two larger lateral references, applied to the lateral quadrants. They can be used in both UPWARD (target above the frame) and DOWNWARD (target below the frame) configurations.


Figure 22

- It also includes two smaller lateral references, intended for the DOWNWARD configuration, when there is limited space between the stereotactic frame and the shoulder, either due to restricted movement or the patient's anatomical characteristics.


Figure 23

9.3.4. MANUAL CALCULATION OF STEREOTACTIC COORDINATES

Stereotactic calculations can be done manually, based on a computerized tomography exam.

 **WARNING:** Find the exam image that best shows the target or region of interest.

Image example (Upward)

Figure 24

- **Determination of the “X” axis**

1. With the help of a Ruler, draw one line connecting points 1 and 3 and another connecting points 7 and 9;

Figure 25

2. Calculate and mark the mid-points (A and B) for each of these lines drawn;


Figure 26

3. Draw a straight line connecting the points (A and B); this line is, therefore, the X axis.


Figure 27

- **Determination of the “Y” axis**

1. With the help of a Ruler, calculate and mark the mid-point (C) of the X Axis (drawn in the previous step);


Figure 28

2. Draw a line connecting points 4 and 6;


Figure 29

3. Calculate and mark the mid-point (D) of the straight line drawn in the previous step;


Figure 30

4. Draw a straight line that passes through the mid-point (D) from the previous step and that passes through the mid-point (C) of the X axis. This line therefore represents the Y axis.


Figure 31

- **Determination of the “Z” coordinate**
- **For “UPWARD” target (positive Z Coordinate)**


Figure 32

- For targets in the RIGHT POSTERIOR quadrant, measure the distance between points 1 and 2 (Z_1);
- For targets in the LEFT POSTERIOR quadrant, measure the distance between points 8 and 9 (Z_3);
- For targets in the RIGHT and LEFT ANTERIOR quadrant, measure the distance between points 5 and 6 (Z_2);
- **For “DOWNWARD” target (negative Z Coordinate)**


Figure 33

- For targets in the RIGHT POSTERIOR quadrant, measure the distance between points 1 and 2 (Z_1);
- For targets in the LEFT POSTERIOR quadrant, measure the distance between points 8 and 9 (Z_3);
- For targets in the RIGHT and LEFT ANTERIOR quadrant, measure the distance between points 4 and 5 (Z_2);

- **Determination of the target's "X" and "Y" Coordinates**

The intersection of the lines drawn previously is the 0 point in "X" and 0 point in "Y." Therefore, based on these lines, we are able to locate the target coordinates, with the following steps:

1. With the help of a Ruler, draw a straight line parallel to the "X" axis, going through the target; draw a straight line parallel to the "Y" axis, going through the target;


Figure 34

2. The intersection of these two new lines should be the target point (point E);


Figure 35

3. Measure the distance (X_1) between point (E) and the "Y" axis. This distance is the "X" coordinate;


Figure 36

4. Measure the distance (Y_1) between point (E) and the "X" axis. This distance is the "Y" coordinate;


Figure 37

NOTE: These coordinates correspond to the virtual coordinates, which should be converted into real coordinates using the method described in the next topic.

• **Conversion of virtual coordinates into real coordinates**

Print images are not always in real scale (1:1), since the properties can be changed to enter a greater number of images in a single print, causing reduction in these images; with this being the case, dimensions found in the print exam (which can be called virtual coordinates) should be processed and then entered into the equipment.

Virtual coordinates should be corrected using a conversion factor, which is found by dividing the Localization's Real Measurement (containing 140mm) by the Localization's Virtual Measurement, found by using a Ruler to measure the distance between points 1 and 3 in the print exam. As shown below:

$$\frac{140}{\text{MVF}} = \text{CONVERSION FACTOR}$$

After finding the conversion factor, the measurements for the virtual coordinates "X", "Y" and "Z" need to be multiplied by the value found, thus providing the real stereotactic coordinates.

• **Example:**

Suppose that the (virtual) measurements found in the CT exam were "X" = 22mm, "Y" = -18mm and "Z" = 18mm. And the measurement found between point 1 and point 3 was 73mm.

$$\frac{140}{73} = 1.92$$

The result (1.92) found is the conversion factor. Multiply the virtual coordinates found by the calculated conversion factor, according to the table below:

Axis	Virtual coordinate	Conversion factor	Stereotactic coordinate
X	22mm	x 1.92	42,24mm
Y	-18mm		-34,56mm
Z	18mm		34,56mm

The location of the target should be shown by the (+) and (-) minus signs:

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Z Coordinate: Top (+) and Bottom (-).

9.3.5. ASSEMBLY OF THE COMPONENTS - COORDINATE TRAYS

⚠ WARNING: Before starting the assembly of the components, it is necessary to set up the Micromar Cranial Headfix fixation system with the frame adapter;

⚠ WARNING: For more information about the Micromar Cranial Headfix fixation system and the frame adapter, refer to the user manual of the associated equipment;

⚠ WARNING: Sterilize the Frame before starting the procedure;

⚠ WARNING: All parts should slide together easily; if this is not the case, see the chapter on lubrication (chapter 13);

⚠ WARNING: If there is any resistance, check that the part is correctly assembled or whether a fixation screw is causing interference;

⚠ WARNING: If the problem persists, contact Micromar or an authorized representative in your region;

⚠ WARNING: Fixation screws for each component should be screwed in before following the step by step instructions below.

- **Identification of the assemblies**

1. Identify which assembly will be mounted first — Assembly A (Coordinate Tray A) or Assembly B (Coordinate Tray B);


Figure 38

⚠ WARNING: The assembly of Set A is identical to that of Set B; the respective trays are identified by having four support feet.

⚠ WARNING: The order in which the assemblies are mounted does not affect the final result.

⚠ WARNING: The parts of Coordinate Tray A must not be mixed with those of Coordinate Tray B, as they have specific assembly adjustments.

- **Ruler Assembly**

2. Place the ruler on the frame;


Figure 39

NOTE: The Ruler can be attached to the four Frame quadrants (right and left sides and anterior and posterior sides), establishing its position according to the target location, as chosen by the neurosurgeon.

3. Ensure that the Ruler is lying entirely against the Frame;


Figure 40

4. Place the two screws and turn clockwise to lock the Ruler in place.


Figure 41

NOTE: When the Upward approach will be used, the fixation screws should be placed in the top part of the Frame, and when the Downward approach will be used, the fixation screws should be placed in the bottom part of the Frame.

- **“Z” Guide Post Assembly**

1. Release all of the fixation screws from the post, so that no interference occurs during assembly;


Figure 42

2. After releasing all of the fixation screws, insert the "Z" Guide Post into the Ruler;


Figure 43

NOTE: When the "Z" coordinate is less than 60mm, the short "Z" Guide Post and the Short "Z" Post can be used.

3. Move the "Z" Guide Post to the desired coordinate;


Figure 44

4. After finding the desired coordinate, use the fixation screws to lock the guide in place, turning it clockwise.


Figure 45

- **“Z” Post Assembly**

1. Release the fixation screw from the “Z” Post, turning it counterclockwise so that there is no interference during assembly;


Figure 46

2. After releasing the fixation screw, insert the “Z” Post (the short “Z” Post can also be used with the short “Z” Guide Post, optional items) into the “Z” Guide Post;


Figure 47

3. Slide the "Z" Post to the desired coordinate;


Figure 48

4. After finding the desired coordinate, use the fixation screws on the "Z" Guide Post to lock the "Z" Post in place, turning it clockwise.


Figure 49

- **Verifying the system's isocenter**

1. Insert the Arc into the Tube;


Figure 50

2. Position the Arc on the desired coordinate and lock it into position by using the fixation screws;


Figure 51

3. Insert the Chariot into the Arc;


Figure 52

4. Position the Chariot on the "80°" angle coordinate and lock it into position by using the fixation screws;


Figure 53

5. On the Chariot's Millimetric Body, go to the "0" coordinate and lock it into position by using the fixation screws;


Figure 54

6. Turn toward "LOCK" and lock in place, so that the Chariot runs no risk of falling;


Figure 55

7. Place the Bushing Guides in the Chariot, screwing them in clockwise;


Figure 56

8. Insert the Cannula in the Bushing Guides attached to the Chariot, until it is against the top Bushing Guide;


Figure 57

9. On the graduated tube, use the sight to check the isocenter according to the cannula model to be used (Figure 58).


Figure 58

• Tube Assembly

1. Remove the Cannula from the Chariot, used to check the isocenter;
2. Check that the fixation screw on the "Z" Post is loose;


Figure 59

3. Insert the Millimetric Tube in the "Z" Post;


Figure 60

4. Position the Tube in the angular position desired, turning it clockwise or counterclockwise;


Figure 61

5. Turn the fixation screw clockwise to lock the Tube into place;


Figure 62

10. SCALES

All of the system's linear scales are graduated in **I** millimeters and all of the system's angular scales are graduated in **II** sexagesimal degrees.


Figure 63

11. TRAYS FOR TRANSPORT / STORAGE / STERILIZATION

For transporting the DoubleArc TM04B Stereotactic System, all parts must be properly placed in their designated locations.

- ⚠ WARNING:** Store in a dry, ventilated and clean environment, protected from directed sunlight and substantial dust;
- ⚠ WARNING:** Three boxes is the maximum amount recommended when piling for transport and storage;
- ⚠ WARNING:** Before using the product, the package should be inspected as to its integrity, in order to guarantee that it is not compromised;
- ⚠ WARNING:** Do not use the material if the package has been damaged and/or violated;
- ⚠ WARNING:** Maximum stacking height, protect from rain, fragile, and this side up are information displayed on the product's packaging and labeling.

12. CLEANING AND STERILIZATION

Cleaning method	
Sanitization	Wet a surgical compress or wipe with enzymatic detergent and apply to the entire surface of the item; rinse using a surgical compress or wipe moistened with water..
Manual Cleaning (immersion)	Immerse the dissembled parts in a container with a solution of cold to warm (up to 40°C) demineralized water and enzymatic detergent (dilute according to manufacturer recommendations). Let the solution act for at least 15 minutes. Using a soft brush, remove residue from the items, paying attention to recesses, screws and lights. Rinse in running demineralized water. Dry using a surgical compress or wipe and clean compressed air.
Automated Cleaning (Ultrasound)	Immerse the dissembled parts in an ultrasonic washing vat with a solution of cold to warm (up to 40°C) demineralized water and enzymatic detergent (dilute according to manufacturer recommendations). Activate the ultrasound for at least 15 minutes. Rinse in running demineralized water. Dry using a surgical compress or wipe and clean compressed air.
Sterilization Method	
VSSP / Autoclave	Sterilization in Saturated Steam under Pressure (SSUP / Autoclave): $134^{\circ}\text{C} \pm 2 \mid \frac{2,14 \text{ kgf/cm}^2}{3 \text{ bar}} \mid 7 \text{ minutes} \mid 100\% \text{ H}_2\text{O} \mid \text{dry for 50 minutes}$
ETO	Sterilization in Ethylene Oxide (ETO): $50^{\circ}\text{C} \pm 10 \mid \frac{-0,3 \text{ bar} \pm 0,1}{-0,3 \text{ kgf/cm}^2 \pm 0,1} \mid \geq 40\% \text{UR} \mid 60 \text{ minutes} \mid 90\% \text{ETO and } 10\% \text{CO}_2 \mid 80 \text{ changes/hour} \mid \geq 48 \text{ hour of aeration.}$

The tables below show the cleaning and sterilization methods for the Stereotactic DoubleArc TM04B :

Item	Cleaning method			Sterilization Method	
	Sanitization	Manual Cleaning (Immersion)	Automatic Cleaning (Ultrasound)	Saturated Steam under Pressure (SSUP / Autoclave)	Ethylene Oxide (ETO)
Fixation Tray	✗	✓	✓	✓	✓
Coordinate Trays A and B	✗	✓	✓	✓	✓
Support for Basic CT	✓	✗	✗	✗	✗
Polyacetal Fix. Screw	✓	✓	✓	✗	✓

Key			
✓	Procedure Permitted	✗	Procedure Prohibited

⚠ WARNING: Items exposed to Creutzfeldt-Jakob disease should be discarded according to local regulations.

⚠ WARNING: If the parameters for cleaning and sterilization are not followed as described here, contamination control may not be guaranteed or damages could be caused to the product.

⚠ WARNING: Do not use abrasive and corrosive agents for cleaning, since they can cause damage to the product.

⚠ WARNING: After cleaning, the components should be immediately dried using a dry compress and compressed air. Pay attention to recesses and screws.

NOTE: Follow the washing machine manufacturer's instructions.

NOTE: Follow the information specified by the detergent manufacturer.

13. LUBRICATION

Before surgery, lubricate the sliding parts and the ends of the fixation screws with sterile silicone or petroleum jelly. Lubricating each part facilitates tightening or loosening the screws and the movement of components during surgery.

14. PREVENTIVE AND CORRECTIVE MAINTENANCE

To guarantee the appropriate functioning and performance of the equipment, it should be returned for preventive maintenance after 18 months.

If there is any defect or malfunction before 18 months have passed, do not use the equipment and immediately contact Micromar or your authorized sales representative.

Micromar is solely authorized to provide technical assistance for the Stereotactic DoubleArc TM04B. After the 18 month period has passed, the equipment must be sent every 14 months to Micromar for assessment to find whether any maintenance or disposal is needed.

PREVENTIVE MAINTENANCE TABLE			
Quality Control	Quality Control	Quality Control	Quality Control
1 st Preventive Maintenance	2 nd Preventive Maintenance	3 rd Preventive Maintenance	4 th Preventive Maintenance

15. SYMBOLS

Symbol	Key
	Not sterile
	Keep out of sunlight
	Manufacturer
	Attention, check documents
	Keep out of rain
	Operating manual, operating instructions
	Date of manufacture
	Serial number
	Stacking limit by number
	Fragile
	This side up
	Lot number

16. WARNINGS AND PRECAUTIONS

- The equipment must only be used with items bearing the same serial number (except for accessories). Replacing components with those of different serial numbers may compromise the equipment's accuracy, causing harm to the patient.
- Accidental movements involving the components of the DoubleArc TM04B Stereotactic System during the procedure may compromise the effectiveness of the operation.
- Avoid drops, twisting, crushing, impacts, or any action that could deform the equipment.
- In any of the above situations, the user must contact Micromar for evaluation.
- Deformation of the equipment structure compromises system accuracy and puts the patient's life at risk.
- Incorrect reading of the equipment scales (due to user error or poor equipment condition) may put the patient's life at risk.
- Ensure all movable parts of the equipment are properly screwed and secured.
- For the procedure, the equipment must be cleaned and sterilized according to the instructions in this manual.
- Immediately after the procedure, the equipment must be cleaned as instructed in this manual.
- The equipment must be sterilized in accordance with the instructions contained in this User Manual; otherwise, there is a risk of bacterial or viral contamination, leading to infectious complications for the patient.
- This product must only be handled by a surgeon specialized in neurosurgery and trained by Micromar.
- The surgeon must consider the patient's physical and pathological conditions when securing the fixation screw of the DoubleArc TM04B Stereotactic System to the patient's skull, in order to avoid injury.
- Failure to follow the instructions contained in this manual and in the user training may result in incorrect use of the DoubleArc TM04B Stereotactic System, putting the patient's life at risk.
- Use Micromar's maintenance service or an authorized representative for proper repair of the DoubleArc TM04B Stereotactic System in case of any deformation.
- For repair or replacement of any component of this equipment, please contact Micromar through customer service or an authorized representative.
- Equipment used to treat patients affected by Creutzfeldt-Jakob disease MUST NOT be reused on other patients. Check local disposal regulations.
- The DoubleArc TM04B Stereotactic System is intended for use only in hospitals and clinics.
- Adverse effects: No secondary or side effects have been reported directly related to the proper use of the product.

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QUALITY CERTIFICATION

This manual is aimed at helping the user understand how to correctly, safely and effectively use the Stereotactic DoubleArc TM04B.. The user must carefully read the user manual that comes with the product before using it.

This product is warranted by Micromar Ind. e Com. Ltda. within the following conditions:

- 1- An 18 month warranty, as from the date of issue of the receipt of sale to the consumer, for replacement of parts or components that may have manufacturing defects, except for those resulting from improper use or installation, done not following the specifications contained in the user manual;
- 2- Warranty pending presentation of this certificate will be offered at no cost to the customer;
- 3- This warranty is null and void in the event of:
 - 3.1- Defect caused by the customer or by third parties other than the manufacturer;
 - 3.2- The product has been changed or repaired by persons or companies not authorized by the manufacturer;
 - 3.3- The recommendations of the user manual are not followed;

Serial Number: _____

Invoice Number: _____

Date: _____

EQUIPMENT INSPECTED BY:

Quality Control

Product	Code	SN	DatE	PREVENT MAINTENANCE TABLE			
Stereotactic DoubleArc TM04B	TM04B	SN: _____	____/____/____	____/____/____ QUALITY CONTROL	____/____/____ QUALITY CONTROL	____/____/____ QUALITY CONTROL	____/____/____ QUALITY CONTROL
				1ST PREVENT MAINTENANCE	2ND PREVENT MAINTENANCE	3RD PREVENT MAINTENANCE	4TH PREVENT MAINTENANCE

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