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
FIXADOR DE CRÂNIO
HEADFIX HF03B

M I C R O M A R

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FIXADOR DE CRÂNIO HEADFIX HF03B

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ENGLISH VERSION

USERS MANUAL

HEADFIX HF03B CRANIAL
FIXATION SYSTEM

PRODUCT DESCRIPTION:

HEADFIX HF03B CRANIAL FIXATION SYSTEM

ONLY USE THIS EQUIPMENT AFTER READING THE INSTRUCTIONS AND RESPECTIVE WARNINGS

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MICROMAR IND. COM. RESERVES THE RIGHT TO MAKE CHANGES TO THIS MANUAL WITHOUT PRIOR NOTICE.
CONTACT OUR SALES DEPARTMENT TO RECEIVE THE LATEST VERSION OF THE USER MANUAL.
THIS MANUAL CONTAINS INFORMATION PROTECTED BY COPYRIGHT.
THIS MANUAL MAY NOT BE REPRODUCED, TRANSLATED OR TRANSMITTED BY ANY ELECTRONIC OR
OTHER MEANS, IN WHOLE OR IN PART, WITHOUT EXPRESS AUTHORIZATION FROM MICROMAR.

ÍNDICE

1. INTRODUCTION.....	33
2. ILLUSTRATIONS.....	33
3. USE OF EQUIPMENT	33
4. RECOMMENDED USE.....	33
5. CRANIAL FIXATION	34
5.1. COMPONENT LIST.....	34
5.1.1. BASE UNIT - HF1000FAA.....	35
5.1.2. ARTICULATED SYSTEM - HF2000FAA.....	36
5.1.3. SKULL CLAMP - HF3000FAA	36
6. EQUIPMENT ASSEMBLY FOR THE SURGICAL PROCEDURE	38
6.1. VERIFYING THE STABILITY OF THE BASE SET PRIOR TO THE SURGICAL PROCEDURE	38
6.2. ASSEMBLY OF THE BASE SET ON THE SURGICAL TABLE.....	39
6.3. ASSEMBLING THE FIXATION PINS IN THE SKULL CLAMP.....	42
6.3.1. PREPARATION TO USE THE FIXATION PIN - ADULT MODEL	42
6.3.2. PREPARATION TO USE THE FIXATION PIN - PEDIATRIC MODEL.....	43
6.4. ADJUSTING THE OPENING DISTANCE OF THE SKULL CLAMP.....	44
6.5. RELEASING THE ADJUSTMENT MOVEMENT OF THE SEMIARC	44
6.6. POSITIONING THE CRANIAL FIXATION SYSTEM ON THE PATIENT'S HEAD	44
6.7. STABILIZATION OF THE PATIENT'S HEAD	47
6.8. INSERTING THE ARTICULATED ARMS INTO THE EQUIPMENT	49
7. DISASSEMBLY OF THE EQUIPMENT AFTER THE SURGICAL PROCEDURE	49
8. CLEANING AND STERILIZATION	51
9. PREVENTIVE AND CORRECTIVE MAINTENANCE.....	52
9.1. ADJUSTMENT OF LEVER PRESSURE	52
10. SYMBOLS.....	55
11. TRANSPORT AND STORAGE	55
12. WARNINGS AND PRECAUTIONS.....	56

1. INTRODUCTION

This manual is aimed at helping the user of the HeadFix HF03B Cranial Fixation System and its accessories to safely and effectively use the equipment.

Before using the equipment (HeadFix HF03B Cranial Fixation System and its accessories), the user should carefully read the manual and gain a total understanding of its contents.

This manual should be kept with the equipment for reference in the event that questions arise.

2. ILLUSTRATIONS

The images in this manual show an illustration of the equipment. There may be some differences in comparison to the physical components. These differences in the illustrations do not interfere with the use of the equipment.

3. USE OF EQUIPMENT

The equipment must be assembled and used as indicated in this manual and for the purpose for which they were created

4. RECOMMENDED USE

The HeadFix HF03B Cranial Fixation System must be used by physicians when the patient's head must be supported and fixed in place for neurosurgical procedures.

The equipment described herein can only be used in combination with other equipment or components that are expressly recognized by Manufacturer as being compatible therewith.

Compatible products:

- Accessories for cranial fixation - HF03B
- Micromar Stereotactic TM03B;
- HeadFix HF02B Cranial Fixation;
- AimNav instruments - Cranial Toolbox - NN00100FAA;
- Disposable Pin for Cranial Fixation (Aluminum);
- Disposable Pin for Cranial Fixation (Titanium);
- Disposable Pin for Cranial Fixation - Stainless steel;
- Reusable Cranial Fixation Pins
- Instrumental for Cranial Fixation.

5. CRANIAL FIXATION

The HeadFix HF03B Cranial Fixation System is equipment that allows for fixation and support of the patient's head as required by the surgeon in neurosurgical procedures.

5.1. COMPONENT LIST

The Cranial Fixation system consists of 3 sets.

ITEM	DESCRIPTION	CODE	QTY.
A	Skull clamp	HF3000FAA	1
B	Articulated sytem	HF2000FAA	1
C	Base unit	HF1000FAA	1


Figure 1

5.1.1. BASE UNIT - HF1000FAA

The Base unit is a set that allows the cranial fixation system to be attached to surgical table.




Figure 2

⚠ WARNING: Before beginning the surgical procedure, ensure that the distances (maximum and minimum) allowed by the Base are compatible with the surgical table to be used during the surgical procedure.

⚠ WARNING: Check the user manual provided by the manufacturer for information on the surgical table.

The Base unit set (HF1000FAA) consists of 3 sets.

ITEM		DESCRIPTION	CODE	QTY.
(A)		Surgical table adaptor	HF1000FBA	1
(B)		Quick Lock lever	HF1000FCA	1
(C)		Stretchers	HF1000FDA	1

5.1.2. ARTICULATED SYSTEM - HF2000FAA

The Articulated System connects the Base unit and the Skull clamp. It allows for 360° pivoting and full flexibility of the equipment in positioning the patient.





Figure 3

5.1.3. SKULL CLAMP - HF3000FAA

The Skull clamp is a set that makes it feasible to fix the patient's head by applying a load distributed across three points around the patient's head.

Minimum opening distance: 80mm


Maximum opening distance: 240mm

Figure 4

⚠ WARNING: Before fix the patient's head verify if space between arms of Skull Clamp with Fixation Pins is larger than patient's head.

⚠ WARNING: Inappropriate adjustment of the cranial fixation system on the patient's head can cause the pins to slip causing damage, such as scalp laceration.

ITEM		DESCRIPTION	CODE	QTY.
A		Fixed arm	HF3100FAA	1
B		Mobile Arm	HF3210FAA	1

6. EQUIPMENT ASSEMBLY FOR THE SURGICAL PROCEDURE

Follow the step-by-step instructions below to assemble the Cranial Fixation system for the surgical procedure:

NOTE: For all assembling step of the Cranial Fixation System, is necessary assistance another person to support the patient's head.

6.1. VERIFYING THE STABILITY OF THE BASE SET PRIOR TO THE SURGICAL PROCEDURE

1. Place the Base unit set on a flat and clean surface **I** and then position the set at a -45° angle in relation to the support surface **II**;


Figure 5

2. Keep the set positioned at -45° in relation to the support surface, then push the lever **a** to lock the entire set;


Figure 6

3. With the entire set locked, apply pressure to the Stretcher set using the palm of your hand **b**, to verify if the pressure of the lever is sufficient to guarantee the set's stability;


Figure 7

⚠ WARNING: During the stability test, the set must not show any movement whatsoever. If any kind of movement occurs, repeat the procedure to adjust the pressure of the lever.

⚠ WARNING: The lever must never be pressed if the Stretcher set is not inserted inside the hole, since it will then damage the working of the Cranial Fixation System.

⚠ WARNING: The pressure of the set comes factory-adjusted.

NOTE: See chapter 9.1 for more information how to adjust the lever pressure.

6.2. ASSEMBLY OF THE BASE SET ON THE SURGICAL TABLE

4. Remove the Skul clamp set from the rest of the equipment, turning the handle on the Articulated System set counterclockwise;


Figure 8

5. First, remove the handle **(c)** from the tube, turning it counterclockwise;


Figure 9

6. After removing the Base unit set's handle, first lift the hex key **(I)** so that it is perpendicular to the handle **(II)**;


Figure 10

7. Insert the hex key into the connection lock screw (left and right) **I**, then turn the handle counterclockwise until the arms are able to be adjusted **II**;


Figure 11

8. Adjust the distance between the centers of the axes according to the distance between the centers of the surgical table connections that will be used for the procedure;


Figure 12

⚠ WARNING: *The Base unit must be positioned on the surgical table in such a way that the set is symmetrical, leaving the tube parallel to the surgical table.*

9. After adjusting the distance between the centers of the axes, lock the arm adjustment mechanism by turning the handle with the hex key clockwise;


Figure 13

10. Return the hex key to its initial position;


Figure 14

11. Return the handle on the Base unit set **I** and then turn it clockwise **II**:


Figure 15

12. Press the lever's safety lock **d** **I** and, with the safety lock still pressed, pull the lever and leave it unlocked **II**:


Figure 16

13. After all of the articulations on the equipment have been released and the Skull clamp has been removed, place the equipment facing the floor, which will make it easier for the user when fixing the patient's head in place;


Figure 17

14. Lock the Base unit set using the fixation resource on the surgical table.

6.3. ASSEMBLING THE FIXATION PINS IN THE SKULL CLAMP

15. First, choose which pin model (Adult **A**) or Pediatric **P**) is most appropriate for the procedure to be performed;


Figure 18

16. After choosing the pin model, they must be sent to the Material and Sterilization Center for sterilization.

⚠ WARNING: See the cranial fixation pins user manual for more informations.

6.3.1. PREPARATION TO USE THE FIXATION PIN - ADULT MODEL

NOTE: For all assembling step of the Cranial Fixation System, is necessary assistance another person to support the patient's head.

17. Insert a Fixation pin with its protective cover into the Dynamometer set;


Figure 19

18. Insert the pins with their protective cover into the fittings that are furthest away from the center of the semiarc, leaving no space between the Fixation pin and the semiarc;


Figure 20

⚠ WARNING: The surgeon must evaluate the patient's conditions, also considering their age, before planning the surgery.

6.3.2. PREPARATION TO USE THE FIXATION PIN - PEDIATRIC MODEL

19. Insert a Fixation pin with its protective cover into the Dynamometer set;


Figure 21

20. Insert the pins with their protective cover into the fittings that are closest to from the center of the semiarc, leaving no space between the fixation Pin and the semiarc;


Figure 22

⚠ WARNING: The surgeon must evaluate the patient's conditions, also considering their age, before planning the surgery.

6.4. ADJUSTING THE OPENING DISTANCE OF THE SKULL CLAMP

21. Remove the protective covers from the Fixation pins;
22. Press the safety trigger **(c)** on the mobile arm **(I)** and, while keeping it pressed the trigger, press the adjustment lock **(d)** and then regulate the distance until the patient's head can be inserted safely **(II)**;


Figure 23

6.5. RELEASING THE ADJUSTMENT MOVEMENT OF THE SEMIARC

23. Turn the semiarc's handle counterclockwise (**(g)**) **(I)**, then adjust the position of the semiarc according to need the surgical procedure **(II)**;


Figure 24

⚠ WARNING: Only turn the semiarc when the lock mechanism is in the unlock (**(g)**), otherwise the internal lock mechanism will be damaged.

6.6. POSITIONING THE CRANIAL FIXATION SYSTEM ON THE PATIENT'S HEAD

NOTE: From this stage forward, another person is needed to provide assistance in supporting the patient's head during the entire procedure of attaching the Skull clamp set to the patient's head.

24. Position the Skull Clamp so that the fixation pins are positioned with the central line of the patient's head;


Figure 25

25. The pins inserted in the semiarc at the time of fixation must be equidistant in relation to the pin inserted into the Dynamometer, resulting in symmetrical and more stable distribution of the load over the patient's head;


Figure 26

26. First, place the Fixed arm over the patient's head (I) and then push the Mobile arm until the Fixation pin is touching the surface of the patient's head (II);


Figure 27

27. Turn the semiarc's handle clockwise (🔒) until the mechanism that locks adjustment movement is in engage;


Figure 28

28. Turn the Mobile arm's handle clockwise (**+**) to apply the load on the patient's head (**I**), after reaching the load required to secure the patient's head, stop rotating the handle (**II**);


Figure 29

The load indicator (**g**) allows the user to view the following load values during the fixation of the patient's head:


Figure 30

N (Newton)	Lbs (Pounds)	Kgf (Quilograma-força)
90 N	~20 Lbs	~9 kgf
180 N	~40 Lbs	~18 kgf
270 N	~60 Lbs	~27 kgf
360 N	~80 Lbs	~36 kgf

Unit conversion table

⚠ WARNING: After actuation the semiarc locking mechanism (**🔒**), it's recommended to check if the locking actually occurred.

⚠ WARNING: The patient's physical and pathological conditions must be taken into consideration prior to applying the load to fix the patient's head.

⚠ WARNING: The inappropriate adjustment of the Skull Clamp on the patient's head could cause the pins to slide, causing serious damage such as laceration of the scalp or falling of the head of the patient, bringing the risk of death.

6.7. STABILIZATION OF THE PATIENT'S HEAD

29. Release the articulation pulling the rod **I** and then turning it clockwise (**II**);


Figure 31

⚠ WARNING: When the articulations are released, ensure that the fixation teeth are separated.

30. Release the turning movement of the articulated System by turning the handle counterclockwise (**II**);


Figure 32

31. Insert the threaded pin **h** into the engage on the Skull Clamp **i** and adjust its position so that the teeth are uniformly connected **I**, then turn the handle clockwise (**II**) to lock attach the Skull clamp to the Articulated system;


Figure 33

⚠ WARNING: Ensure that the correct placement of the fixation teeth during fixation.

32. Pull the rod **f** **I** then turning it counterclockwise (**II**), locking the articulation of the Articulated system **II**;


Figure 34

33. Adjust the position according to need of the surgical procedure **I** then push the lever until the locking mechanism is engage **II**;


Figure 35

⚠ WARNING: The lever must never be pressed if the Stretcher is not inserted inside the hole, since it will then damage the working of the Cranial Fixation System.

⚠ WARNING: The pressure of the Quick lock lever comes factory-adjusted.

⚠ WARNING: Make sure that at the time of fixing the locking teeth are correctly engaged.

⚠ WARNING: After activating the semiarc's locking mechanism (**II**), it's recommended to check if the lock actually happened.

⚠ WARNING: The patient's physical and pathological conditions must be taken into consideration prior to applying the load to fix the patient's head.

⚠ WARNING: The inappropriate adjustment of the Skull Clamp on the patient's head could cause the pins to slide, causing serious damage such as laceration of the scalp or falling of the head of the patient, bringing the risk of death.

6.8. INSERTING THE ARTICULATED ARMS INTO THE EQUIPMENT

34. Insert the threaded pin **(j)** into the connection on the Skull clamp **(k)** and adjust its position so that the teeth are uniformly connected **(I)**, then turn the Articulated arm's handle clockwise to attach it to the Skull clamp set **(II)**.


Figure 36

⚠ WARNING: Ensure that the correct placement of the fixation teeth during fixation.

7. DISASSEMBLY OF THE EQUIPMENT AFTER THE SURGICAL PROCEDURE

Follow the step-by-step instructions below to disassemble the Cranial Fixation system after the surgical procedure:

NOTE: For all disassembling step of the Cranial Fixation System, is necessary assistance another person to support the patient's head.

1. Turn the Articulated system's handle counterclockwise (**(II)**), to separate the Skull clamp set from the rest of the equipment;


Figura 37

2. Press the lever's safety lock **(I)** and, while keeping the safety lock pressed, pull the lever and leave it unlocked **(II)**:


Figura 38

- Place the equipment without the Cranial Fixation set downward facing;


Figura 39

- Turn the Mobile arm's handle counterclockwise (●) (I) to reduce the load applied to the patient's head (II);



Figura 40

- Press the safety trigger on the Mobile arm (I) and, while keeping it activated, press the adjustment lock and then regulate the distance until the patient's head can safely be removed (II);


Figura 41

- Send the entire equipment for cleaning as specified in this manual.

⚠ WARNING: See chapter 8 for more information on the correct cleaning process.

8. CLEANING AND STERILIZATION

The table below shows the methods of cleaning and sterilizing the equipment:

Cleaning Method	
Sanitation	Wet a surgical compress or wiper with enzymatic detergent and apply over the entire surface of the item; rinse with surgical compress or wiper moistened with water.
Manual Cleaning (immersion)	Immerse the disassembled parts in a container with cold to warm demineralized water (up to 40°C) with enzymatic detergent (dilution as recommended by the manufacturer). Allow the solution to act for at least 15 minutes. With a soft brush remove the residue from the items, taking care of the recesses, threads and lumens. Rinse with running demineralized water. Dry with surgical compress or wipe and clean compressed air.
Cleaning Automated (Ultrasound)	Immerse the disassembled parts in an ultrasonic washer tank with cold to warm demineralized water solution (up to 40°C) with enzymatic detergent (dilution as recommended by the manufacturer). Operate the ultrasound for at least 15 minutes. Rinse with running demineralized water. Dry with surgical compress or wiper and clean compressed air.
Sterilization Method	
Steam	Sterilization in Saturated Vapor under Pressure (Steam): 134°C ±2 $\frac{2,14 \text{ kgf/cm}^2}{3 \text{ bar}}$ 7 minutes 100% H ₂ O 50 minutes of drying
ETO	Sterilization in Ethylene Oxide (ETO): 50°C ±10 $\frac{-0,3 \text{ bar} \pm 0,1}{-0,3 \text{ kgf/cm}^2 \pm 0,1}$ ≥40%UR 60 minutes 90%ETO e 10%CO ₂ 80 exchanges/hour ≥ 48 hours of aeration

Item	Code	Cleaning Method			Sterilization Method	
		Sanitation	Manual Cleaning (immersion)	Cleaning Automated (Ultrasound)	Saturated Steam Under Pressure (Steam)	Ethylene Oxide (ETO)
Headfix HFO3B Cranial Fixation System	HFO3B	⊗	⊙	⊙	⊗	⊗
Subtitle		⊙	Allowed Procedure		⊗	Prohibited Procedure

- ⚠ **WARNING:** Items exposed to Creutzfeldt-Jakob disease should be discarded according to local regulations;
- ⚠ **WARNING:** If the parameters for cleaning and sterilization are not followed as described herein the contamination control may not be guaranteed or damage may be caused to the product;
- ⚠ **WARNING:** Do not use abrasive and corrosive agents for cleaning, since they can cause damage to the product;
- ⚠ **WARNING:** After cleaning, the components should be immediately dried using a dry compress and compressed air. Pay attention to recesses and screws;

- NOTE:** Follow the manufacturer instructions for the washing machine;
- NOTE:** Follow the information specified by the detergent manufacturer.

9. PREVENTIVE AND CORRECTIVE MAINTENANCE

To guarantee the appropriate functioning and performance of the equipment, it must return to preventive maintenance **after 18 months**.

If there is any defect or malfunction before 18 months have passed, do not use the equipment and immediately contact Micromar or your authorized sales representative.

Micromar is the only authorized technical service provider for the Headfix HF03B Cranial Fixation System.

After the expiration of the 18-month period, it's need sent the equipment to Micromar, and the **every 14 months**, for assessment and to perform any maintenance or disposal.

PREVENTIVE MAINTENANCE TABLE			
			
Quality control	Quality control	Quality control	Quality control
1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

9.1. ADJUSTMENT OF LEVER PRESSURE

 **WARNING:** Check the lever tension at regular intervals.

 **WARNING:** Adjustment of lever pressure must only be done when there is a real need to do so, since lever pressure already comes factory-adjusted.

1. Press the lever's safety   and, while keeping the safety lock pressed, lock pull the lever and leave it unlocked .


Figure 42

2. Remove the handle from the tube, turning it counterclockwise  and then remove the Base unit .


Figure 43

3. After removing the Base set's handle, first lift the hex key **I** until it's perpendicular to the handle **II**;


Figure 44

4. Insert the hex key into the hexagonal profile of the Quick lock lever **I** after turning the handle counterclockwise to release the pressure adjustment movement of the lever **II**;


Figure 45

5. Hold the lever all the way up **I** and then turn the handle clockwise manually to apply more pressure (**+**) or counterclockwise to apply less pressure (**-**) **II**;


Figure 46

⚠ WARNING: After adjusting the lever pressure, ensure that the lever's movement is smooth yet firm.

6. After adjusting the lever's pressure, press the lever until it locks **I**, then insert the hex key into the hexagonal profile of the Quick lock lever set and turn the handle clockwise to lock the lever's pressure adjustment movement **II**;


Figure 47

- ⚠ WARNING:** The lever must never be removed from the set;
- ⚠ WARNING:** The lever must never be pressed if the Stretcher set is not inserted inside the hole, since it will then damage the working of the Cranial Fixation System.
- ⚠ WARNING:** Ensure that the backrest profile of the lever is correctly placed on the support pins;
- ⚠ WARNING:** If the lever is not connected to the support pins, carefully reinsert it, ensuring that the backrest profile of the lever is correctly placed on the support pins.
- ⚠ WARNING:** Excessive pressure on the lever could damage the locking mechanism.

10. SYMBOLS

Symbols	Subtitle
	Keep away from sunlight/Keep away from heat
	Manufacturer
	Attention: Consult the instructions for use
	Keep dry/Keep away from humidity
	Authorized representative in the European Community
	Consult the instructions for use
	Date of manufacture
	Serial number
	Stacking limit
	Fragile
	This way up
	Lot number

11. TRANSPORT AND STORAGE

The product must be transported and stored in the packaging provided by Micromar.

12. WARNINGS AND PRECAUTIONS

- Accidental movements involving cranial fixation system components during the procedure can compromise the efficacy of the procedure.
- Avoid dropping, twisting, crushing, impacting or any action that will deform the equipment.
- In the event that one of the above occurs, the user should contact Micromar for assessment of the equipment.
- Incorrect reading of equipment scales (user error in reading or because of poor equipment condition) can put the patient's life at risk.
- Screw in and correctly attach the movable parts of the equipment.
- To perform the procedure, the equipment must be cleaned according to the instructions in this manual.
- Immediately after ending the procedure, the equipment must be cleaned according to the instructions in this manual.
- This product is solely made for use in surgical centers.
- This product may only be manipulated by a surgeon specialized in neurosurgery with experience performing cranial fixation.
- The surgeon must consider the patient's physical and pathological conditions when attaching the pin for cranial fixation along with the equipment (HF03B) to the patient's head, in order to prevent damage to same.
- If the instructions contained in the manual user supplied by Micromar are not followed, this could result in incorrect use of the Cranial fixation system and its accessories, placing the patient's life at risk.
- Use the maintenance service provided by Micromar or an authorized representative to correctly repair the system when there are any problems with its functioning.
- To repair or replace any component of this equipment, you must contact Micromar through customer service or through any authorized representative.
- Equipment used to treat patients suffering from Creutzfeldt-Jakob disease **MUST NOT BE** reused for other patients.
- Disposal of the product should be done according to local laws.
- The basic and optional composition have a useful life of 5 years, as from the date of issuance of the invoice.

Avenida Marginal ao Córrego da Serraria, nº168 • SP • Brazil • Phone: (+55 11) 4057-5722

QUALITY CERTIFICATION

CRANIAL FIXATION SYSTEM AND ACCESSORIES

This product is guaranteed by Micromar Ind. e Com. Ltda. under the following conditions:

- 1- An 18 month warranty, as from the date of issue of the receipt of sale to the consumer, for replacement of parts or components that may have manufacturing defects, except for those resulting from improper use or installation, done not following the specifications contained in the user manual;
- 2- Warranty pending presentation of this certificate will be offered at no cost to the customer;
- 3- This warranty will be null in the event:
 - 3.1- Of defect caused by the customer or by third parties other than the manufacturer;
 - 3.2- The product has been changed or repaired by persons or companies not authorized by the manufacturer;
 - 3.3- The recommendations of the user manual are not followed;
- 4- EAny defects resulting from natural wear to the product or caused by agents of nature are not covered by this warranty. After expiration of the stipulated term, the HeadFix and Accessories should be sent to Micromar, for assessment of any maintenance and/or disposal of the product.
- 5- After the warranty has expiry and the every 14 months, the Cranial fixation system and the accessories must be sent to Micromar for evaluation of the equipment to detect possible maintenance or disposal of the equipment.

*** This manual and certificate should be kept with the product and they may be sent along with the equipment each time it is sent for preventive or corrective maintenance.**

Equipment inspected by:

Quality Control

Product	Code	SN	Date	PREVENTIVE MAINTENANCE TABLE			
HeadFix HF03B Cranial Fixation System	HF03B	SN.: _____	____/____/____	____/____/____	____/____/____	____/____/____	____/____/____
				Quality control	Quality control	Quality control	Quality control
				1st Preventive maintenance	2st Preventive maintenance	3st Preventive maintenance	4st Preventive maintenance

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