





GETTING STARTED

Welcome to Inner Armor - Core!

1. Download the App

- Visit App store  or Google Play store  and search for "Inner Armor"



2. Log-in using your credentials

- Locate your welcome email and search for your password to use along with your email address

A screenshot of the Inner Armor login form. The form has a dark background with white text and input fields. At the top is the "INNER ARMOR" logo. Below it is an "Email" field containing "example@samplemail.com". Underneath is a "Password" field with a toggle icon for visibility. To the right of the password field is a "Forgot Password?" link. At the bottom left is a "Stay logged in?" checkbox. At the bottom right is an orange "LOGIN" button.

3. Select/Click Inner Armor CORE




4. Get breathing belt

- Hold power button on breathing belt for 3 -5 seconds to turn on. Green light should appear to indicate power
- Put breathing belt on so that the sensor is directly in front of you with the belt fastened and tightened to your comfort
- Sensor should rest directly above your belly button
- The Inner Armor helmet should be facing forward so you can see the green light




5. Press start session

- Make sure you are in a comfortable seated position with good posture for the first session




6. Connect the breathing belt via Bluetooth

- Back in the app select Synchroni Breathing Belt





- Search for the corresponding breathing belt number to appear under “Available Devices”
 - The last four digits on the top of the belt should correspond the four letters/numbers on the screen
- Select the number to pair the breathing belt (note in android you must pair the belt within the device settings, before paring in the app)





7. Once the belt is connected the warmup screen automatically appears with a breathing pacer of 10 seconds
- At the end of the 10 seconds the setting session will begin automatically

Please watch the next video for any overview of the session.

- • Check out our FAQs at www.forgeinnerarmor.com/core