



CMP[®] T Cell Medium

Chemically Defined, Serum-Free, and Xeno-Free Medium cGMP-Compliant Manufacture



- T cells direct from PBMCs — no beads, no sorter
- Serum-free, Chemical-defined, cGMP-compliant manufacture
- Day-14 expansion >500x with high naïve and Central Memory T-Cell
- Supports T cell growth in multiple expansion platform
- Custom sizes and packaging available upon request



- A serum-free T-cell expansion system designed for direct PBMC culture, delivering robust proliferation comparable to serum-containing competitors.



- Promotes CD8 enrichment and enhances Naïve/TCM phenotypes for superior T-cell expansion.



Product	Size
T-Cell Expansion Medium	Bottle : 500 mL / 1,000mL
	Bag : 250 mL / 500mL / 1,000mL
T-Cell Expansion Medium-Supplement A	1 vial
T-Cell Expansion Medium-Supplement B	1 vial